



Immedia Rescue

EMERGENCY EVACUATION SLIDE SYSTEM






Immedia RescueSheet is an evacuation sheet that remains in place under the mattress. In case of emergency, it can be used to quickly transfer a bedridden patient including the entire mattress out of the room and down the stairs. The material is waterproof.



Rescue emergency evacuation slide sheet.

INTENDED USE

- Emergency evacuations
- Transfer out of bed, and slide out to safety

DESCRIPTION	WIDTH X LENGTH		SPECIFICATIONS
Immedia RescueSheet	CM	INCH	MATERIAL: Nylon STRAP: Polypropylene BUCKLE: Acetal
~ Standard	75 x 200	29.5 x 78.7	
~ Extra wide	85 x 200	33.1 x 78.8	
WASHING INSTRUCTIONS			
AUTOCLAVE: 130 °C. Can be wiped off with soapy water or disinfectant Do not use fabric softener or bleach.			
		  	 

Immedia HipTurner

SPECIAL AID FOR TURNING AND STATIC SUPPORT IN BED

Immedia HipTurner is especially designed for turning a patient and supporting them in a static position on their side, while also creating space for placing a bed pan under the patient, for personal hygiene, ulcer care or dressing. HipTurner makes it possible for the caregiver to maintain a good working position thus minimizing strain.








Turning with help of Immedia HipTurner can be performed manually or in combination with a hoist.



Immedia HipTurner supports turning and static positioning of the patient on one side.

INTENDED USE

- Turning in bed
- Support static positioning on one side
- Facilitate placing bed pan, personal hygiene, ulcer care or dressing
- Can be used in combination with hoist

DESCRIPTION	WIDTH X LENGTH X HEIGHT		SPECIFICATIONS
Immedia ~ HipTurner (sling 25 x 80 cm)	CM 40 x 20 x 4,5	INCH 15.7 x 7.9 x 1.8	MATERIAL: Nylon STRAP: Polyester COVER: Nylon with PU-coating FOAM: Polyethylene
WASHING INSTRUCTIONS Do not use fabric softener or bleach. Remove the foam pad before cleaning.			
<div></div>			<div><div>max 200 kg</div><div></div><div>max 440 lbs</div></div>