



## ROHO® CONTOUR SELECT® Cushion

Designed for positioning and posture, the CONTOUR SELECT cushion keeps you centered in the middle of the cushion, with your pelvis back in the wheelchair. The ISOFLOR® Memory Control locks air into four separate compartments, providing stability and positioning by minimizing side-to-side and front-to-back motion.

- ✓ **Pressure injury prevention and treatment**
- ✓ **In-place individualization with ISOFLOR® Memory Control**
- ✓ **Four independent compartments for customized fit**
- ✓ **Stability**
- ✓ **Positioning**

### Dimensions

Width: 14.5 in. to 24.25 in. / 36.5 cm to 62.5 cm

Height: Approx. 2.25 in. to 4.25 in. / 5.5 cm to 10.5 cm

Range due to varying heights of cells

Weight: Approx. 3.0 lbs. / 1.4 kg

### Weight limit

No weight limit. However, the cushion must be properly sized to the individual.

### Limited warranty

Cushion: 36 months  
Cover: 6 months



### Our evidence shows the CONTOUR SELECT is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Need added positioning and stability
- Need leg alignment and pelvis stabilization
- Need additional relief under ischial tuberosities (ITs)