

 QUICKIE®


QXi

Innovation, durability, and adjustability
in an ultra-lightweight value package



Live Without Limits

The Quickie QXi is a durable, ultra lightweight folding wheelchair with quick and intuitive adjustments at an outstanding value.



Adjustable Axle



Low Profile, Single Post Armrest



Three-Arm Crossbrace

Quick and Intuitive Adjustments

Most adjustments can be made in less than one minute, using only one tool.

- Center of Gravity
- Seat to Floor Height
- Camber
- Axle Spacing
- Back Height

Simple Option Mix

- Options include the most frequently used armrests, wheels, tires, and seat widths and depths
- New Swing-in/Swing-out Hangers with simple footrest adjustments and intuitive release levers
- Transit option

Durable Frame Design

- Ultra lightweight and durable aluminum frame with a 265 lb. weight capacity
- Three arm cross brace that sits inside frame for a more rigid ride



SPECIFICATIONS Quickie QXi

HCPCS Code	K0005
Frame Style	Folding – Traditional
Transport Weight	17.6 lbs
Overall Weight	28.8 lbs
User Weight Capacity	265 lbs
Seat Width	14" to 20"
Seat Depth	14" to 20"
Back Angle Adjustment	-15° to 3°
Front Seat-to-Floor Height	12.5" to 19.5"
Rear Seat-to-Floor Height	12" to 19"
Hanger Angle	60°, 70°, 80° & Elevating
Warranty	Limited Lifetime