


The new Invacare® Premier Series slings balances solutions for any transfer. We say that with confidence because these slings were developed in collaboration with Occupational and Physical therapists. Invacare Premier slings are available in a wide variety of styles and can be adjusted to ensure a comfortable and secure fit for the resident. Clearly labeled slings and easy attachment points make it simple for the caregiver to quickly attach the sling and lift the resident.

**Application:** For use with Invacare® Reliant 350, 450 and 600 series power lifts, the Invacare® Premier Series Jasmine™ and Roze™ lifts.

**Models:** ROZE, JASMINE, RHL450-1, RPS350-1, RPS350-2, RPL450-1, RPL450-2, RPL600-1, RPL600-2



## UNIVERSAL HIGH SLING PLUS



The Universal High Sling supports the whole body and head in transfers to/from seated or lying position. It offers a comfortable recline for those with limited head, trunk or hip control and allows arms to hang outside the sling.


*Compatible with the following lifts*

- Jasmine™
- Reliant™ 450, 600

### PRODUCT NUMBER/MATERIAL

SIZE	A	B	C	D	POLYESTER FABRIC
S	54.4"	33.1"	34.7"	13.4"	2485960
M	59.5"	37.4"	36.6"	17.7"	2485961
L	62.3"	43.3"	38.2"	22.1"	2485962
XL	62.6"	46.9"	38.2"	24.8"	2485963

## STAND ASSIST SLING



The Stand Assist Sling offers a comfortable and safe transfer for those who have some weight bearing ability, but cannot stand independently. The sling is well padded and features non-slip material on the side that touches the resident to assure a comfortable transfer position.


*Compatible with the following lifts*

- ROZE™
- Reliant™ 350

### PRODUCT NUMBER/MATERIAL

SIZE	A	B	C	POLYESTER FABRIC
S	10.2"	33.1"	34.7"	2484685
M	13.0"	35.3"	46.1"	2484686
L	14.2"	37.4"	50.0"	2484687
XL	14.2"	40.2"	71.3"	2484688

## TRANSFER STAND ASSIST SLING



Also designed for optimal use with our Roze and Reliant 350 lifts, the Transfer Stand Assist Sling is a high-quality sling which provides simple seated transfer for residents with minimal weight bearing ability.

*Compatible with the following lifts*

- ROZE™
- Reliant™ 350

### PRODUCT NUMBER/MATERIAL

SIZE	A	B	C	POLYESTER FABRIC
S	33.9"	33.1"	9.5"	2484689
M	40.2"	35.3"	10.2"	2484690
L	46.5"	37.4"	13.0"	2484691
XL	48.5"	40.2"	13.8"	2485117