
Akron Streamline Changing Table

Instruction for Use


ARJOHUNTLEIGH



GETINGE GROUP



7. Technical Data

General	
Safe working load	225kg
Product weight (approx.)	83kg approx.
Audible noise	54dB
Operating conditions	
Temperature	5°C to 40°C
Relative humidity	20% to 90% at 30°C, non-condensing
Atmospheric pressure	700hPa to 1060hPa
Upholstery material	Mundial Vinyl: 88% PVC - 12% Cotton Diabolo Vinyl: 84% PVC - 16% Cotton

Electrical data		
Power input	230V ~ 2A 50/60Hz	120V ~ 4A 60Hz
Duty cycle	10% (2 min. on, 18 min. off)	10% (2 min. on, 18 min. off)
Electric shock protection	Class II Type B	Class II Type B
Liquid ingress protection	IPX4	IPX4
EMC	Complies with EN 60601-1-2	
	Classified by Underwriters Laboratories Inc. with respect to electric shock, fire, mechanical hazards and other specified hazards only in accordance with the List of Standards and Certificate.	

Environmental Protection	
	<p>This product is marked with the crossed out wheeled bin symbol to indicate that it is covered by EU Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). This symbol indicates that all electrical and electronic products must be taken to a dedicated collection point at end of their working life. Contact ArjoHuntleigh or your local authority for advice on proper disposal.</p>
	<p>Incorrect disposal of this equipment and its component parts, particularly gas springs, actuators and other electrical devices, may produce substances that are hazardous to the environment. To minimise these hazards, contact ArjoHuntleigh for information on correct disposal.</p>