


How to buy

1. **Measure** your client to determine the size you need. (See dimension chart for help.)
2. **Select** the appropriate Mobile Stander and color.
3. **Choose** the accessories you need.



Retail Order Form
August 2017

All Rifton products are made to order and are non-returnable and non-exchangeable.

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.				
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25x25	25x25½	32x29	39x33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without big wheels	19½	21	24	28½
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lbs)	28	30	40	65
Max. working load (lbs)	50	50	120	175



The **conversion kit** converts our smallest stander, the K110, to the next size up, the K130, for less than one quarter the cost of a new stander. Request from customer service.

Conversion kit K116



		K110 mini	K130 small	K150 medium	K170 large
Mobile Stander frame includes body support with padding, seat pad, straps and fleece. Specify color: raspberry, blue or lime		K101	K103	K105	K107
		RAS <input type="checkbox"/>	RAS <input type="checkbox"/>	RAS <input type="checkbox"/>	RAS <input type="checkbox"/>
		BLU <input type="checkbox"/>	BLU <input type="checkbox"/>	BLU <input type="checkbox"/>	BLU <input type="checkbox"/>
		LIM <input type="checkbox"/>	LIM <input type="checkbox"/>	LIM <input type="checkbox"/>	LIM <input type="checkbox"/>
Optional accessories	Wheels (pair) (includes brakes)	K113	K113	K115	K117
	Mounting bracket	K133 <input type="checkbox"/>	K133	K173	K173
	Arm prompts (pair) (require mounting bracket)	K512	K512	K512	K512
	Hand loops (require mounting bracket)	Single	K578	K578	K578
		Pair	K579	K579	K579
	Communication tray (requires mounting bracket)	K518	K518	K518	K518
	Additional seat pad	Small	N/A	K135	K135
		Large	N/A	K155	K155
	Sandals with wedges (pair)	Small	K114	N/A	N/A
		Medium	N/A	K154	N/A
		Large	N/A	N/A	K174