# THE EVOLUTION ARM TROUGH SYSTEM.

The Evolution Arm Trough System has been designed for the person that cannot use a traditional walker handle. It provides a more stable platform, allowing users to remain upright.

Evolution's innovative Arm Trough System is available only for Evolution's four Evo models.

## FEATURES.

Provides a more stable platform; allowing user to remain upright.

5 points of adjustability for custom fitting.

No extra tools required; fully adjustable using built-in tools.

Can be used as a pair or single side with a regular handle.

Can be adjusted to fit hemiplegics.

Otto Bock arm pads can be interchanged or combined.

Removable 7" midsection for extra low settings.

Maximum total weight capacity with two Arm Trough 400lbs.

Arm Trough System available only for Evo Regular, Evo Low, Evo Wide and Evo Low.

Optional flexible back strap paid available.

## SPECIFICATIONS.



### Low & Low wide.

24" Width (Low)

28" Width (Low Wide)

41" Height (with midsection)

7" Removable Midsection

400 lbs Weight Capacity

### Regular & wide.

24" Width (Regular)

28" Width (Wide)

34" - 48" Elbow Pad Height Range 34" - 50" Elbow Pad Height Range

41" Height (with midsection)

7" Removable Midsection

400 lbs Weight Capacity



Port Coquitlam, BC V3C2J9

**EVOLUTION TECHNOLOGIES USA INC.** 2657 Delta Ring Road Ferndale WA 98248