

Pediatric Walkers

- Pediatric walking aids feature “grow-with-you” height adjustability
- Models with wheels maneuver well over normal floor obstacles and are ideal for outdoor play areas
- Guardian pediatric walkers have a design to meet everyone’s needs, even the smallest child



G07749



G07751



G07782

Children’s Folding Walker

Lightweight aluminum frame maintains “A” frame design to provide greater balance. Children’s folding walkers have comfortable vinyl hand grips designed for a child’s hand.

Children’s Non-Folding Walker

Lightweight, aluminum frame with proportionally large width to length for a wide stance and greater stability. Large rubber tips come standard to provide long braking action as this model comes with 3" (7.62 cm) front wheels.

Strider Prescription Walker

The infinite configurations accommodate the special needs of children with disabilities. Handles adjust without tools to allow an infinite number of positions for desired posture. Strider walkers have an extra wide base for maximum stability, safety and balance.

Item Number	Description	Wheels Inch (cm)	Weight Capacity Lbs (kg)	Approximate User Height Inch (cm)	Width Inside Hand Grip Inch (cm)	Height Adjustment Inch (cm)	Folded Width Inch (cm)	Pkg
Strider Prescription Walker								
G07782	Folding, Red	5 (13)	90 (41)	2'2"-3'6" (66-107)	Adjustable	14-19.5 (36-50)	9.5 (24)	Each
G07783	Folding, Green	5 (13)	175 (79)	3'-5'2" (91-157)	Adjustable	24-36 (61-91)	9.5 (24)	Each
Children's Folding Walker								
G07749	Children's	None	200 (91)	4'-4'8" (122-142)	12.75 (32)	24.5-28.5 (62-72)	8 (20)	Each
G07750	Toddler	None	175 (79)	3'4"-4' (102-122)	10.75 (27)	20.5-24.5 (52-62)	8 (20)	Each
Children's Non-Folding Walker								
G07751	Tyke	3 (8)	80 (36)	2'8"-3'6" (81-107)	13 (33)	15.5-20.5 (39-52)	N/A	Each