

Accessories



Sanitary Collection Pan with Lid

Model no. 10230

Available for all Ocean models

Lid only - model no. 1541289




Hygiene Back Support

Contoured back support is easy to use and clean, as well as ergonomically shaped to fit all users.

Model no. 1575299

Available for Ocean Standard, SP, VIP, E-VIP only



Special Soft Seat

Soft, ergonomic watertight surface with a non-slip textured surface and Gore-Tex® membrane.



Hygiene Recess

Soft Seat

With hygiene recess.

Model no. 16341

Available with Ocean Standard, SP, XL, E-VIP only

Model no. 1470075

Available for Ocean VIP only



Toileting Recess

With toileting recess.

Model no. 16342

Available for Ocean Standard, SP, XL, E-VIP only



Insert for Soft Seat

For soft seat with hygiene recess.

Model no. 16343

Available for all Ocean models



Calf Support with Foot Support

Left or right side. Adjustable in height and angle. Upholstered pads.

Model no. 1546705 (Left), 1546439 (Right)

Available for all Ocean models



Ocean Laterals

Adjustable lateral supports to provide trunk support where upper body control is limited or in question.

Model no. 1535077

Available for all Ocean models



1470078

Safety Bar

For additional security. Easy assembly. Padded grip.

Model No. 10167-10

Available for Ocean Standard, SP, VIP, E-VIP only

Model No. 1470078

Available for Ocean XL only