



Groupe
Medicus

24"/61 cm, knurled stainless steel grab bar. Features HealthCraft's innovative Easy Mount 9-hole flange design for greater mounting flexibility.

Easy Mount™ Flange - Unlike 3-hole flanges that limit screw placement, the innovative Easy Mount 9-hole flange allows the grab bar to be installed at varying degrees.

Compatible Flange Design - Also compatible with specialized grab bar fasteners and anchors that mount directly behind drywall, tile and other surfaces should stud spacing not permit direct mounting.

Will Never Rust - Stainless steel construction ensures that these grab bars will never rust. A simple wipe with a damp cloth is all that is required to maintain a clean attractive finish.

Concealed Screws - Flange covers conceal the screws for an attractive finish.

ADA Compliant - All HealthCraft™ grab bars meet or exceed the guidelines of the Americans with Disabilities Act (ADA).

Material:	Knurled stainless steel.
Diameter:	1.25"/32mm
Overall Length:	24"/61cm
Clearance:	Distance between the rail and the wall is 1.5"/38mm
Product Weight Capacity:	400 lbs. / 181 kg.
Warranty:	One year limited.