



















Rear Glide Brakes with 3" Wheels

Rear Glide Brakes

Item #	UOM
10106	1/pr

- For use on back of wheeled walker.
- 3" wheel glides along surface until weight is applied to hand grip, which activates rubber-tip brake.
- Allows for 8 height adjustments.
- 350 lb. Weight Limit.

Item #	UOM
10104	1/pr

- For use on back of wheeled walker.
- Nylon tip glides along surface until weight is applied to hand grip, which activates rubber-tip brake.
- Allows for 8 height adjustments.
- 350 lb. Weight Limit.

specifications	
Weight Limit	350 lbs

specifications	
Weight Limit	350 lbs.