



FEATURES

- Contoured seat prevents child from sliding forward and improves seating posture
- 95° seat to back angle
- Double-strap trunk support and a positioning belt
- Two sets of straps secure Seat2Go to any chair. When not in use, the seat straps can be stored under the seat behind the clear flap and the back straps in the back pocket
- Can be used as an insert for a single or double stroller
- Weighs only 5 lb



Seat2Go with optional Height Adjustable Headrest and Depth Adjustable Abductor

ACCESSORIES



Height Adjustable Headrest



Depth Adjustable Abductor

ORDER FORM **FAX TO 905-787-8550 OR SUBMIT ONLINE**

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

| BASE UNIT | CRS 2000 | CRS 3000 |
|----------------------------|----------|----------|
| Product Size | Small | Medium |
| Product Width | 8" | 12" |
| Product Depth | 9" | 12" |
| Product Height of Back | 12" | 14" |
| Product Weight | 5 lbs | 7 lbs |
| ACCESSORIES | | |
| Height Adjustable Headrest | CRS 8000 | CRS 8000 |
| Depth Adjustable Abductor | CRS 8025 | CRS 8025 |