



Etac Swift Freestanding toilet seat raiser

**Back support**

84005069

Material: Polypropylene

**Soft Back**

84005065

Material: Polyethylene

**Soap dish/ Lid holder**

80209266

Material: Polypropylene

**Pan**

without lid 84005066

with lid 84005068

Material: Polypropylene

**Soft Seat**

84005064

Material: Polyethylene

Accessories - Etac Kaskad Freestanding toilet seat

**Bucket with handle and lid**

80302026

Material: Polypropylene

**Bucket holder**

80302025

For bucket with handle.

Material: powder coated steel

**Plastic Bucket**

80302023

Fits without bucket holder.

Material: Polypropylene

Accessories - Etac Cloo toilet seat raiser

**Arm supports**

Pair: 80301222.

Cloo without arm support can be supplemented with arm supports.

Assembly without tools

Material: Polyamide

**Heightening spacers**

4 pcs: 80301221

To obtain 14 cm (5½") seat height, a set of extra spacers is required.