

Please mark clearly with an X the items you wish to purchase.  
Email form to NuVision Rehab Group at [sales@nuvisionrehab.com](mailto:sales@nuvisionrehab.com) or fax form to 403-253-1744.  
Contact our customer service team at (877) 252-6655.

|                  |                      |
|------------------|----------------------|
| Account #: _____ | Ship to: _____       |
| Date: _____      | Name: _____          |
| PO #: _____      | Address: _____       |
| Buyer: _____     | City: _____          |
| Attention: _____ | Phone: (     ) _____ |

## EIO Push Chair Features:

- Removable hip lateral and trunk lateral supports
- Height adjustable padded headrest
- Back rest angle adjusts from 25° to 75°
- Footwell and angle adjustable push handle
- Adjustable 5-Point padded chest harness and shoulder pads as well as a padded 3-point locking pelvic belt
- Removable sun canopy with window and pocket
- Solid filled 9" front and 11" rear wheels require no air
- Swivel front wheels with lock-out feature add security on rough terrain
- Angle adjustable foot well with a 2-position foot-block.
- Rear spring suspension for comfortable ride
- Easy access one-touch parking brake
- Lightweight aluminum frame
- Shopping basket included



EIO Push Chair



Liners/Sitters

Includes  
2-Position  
Foot-Block



## Optional Accessories:



Rain Cover



For kids who need mild to moderate positioning, the Size 1 and 2 Soft-Touch® Liners and Size 1 and 2 Soft-Touch® Sitters attach to the stroller.

Please mark clearly with an X the items you wish to purchase.  
Email form to NuVision Rehab Group at [sales@nuvisionrehab.com](mailto:sales@nuvisionrehab.com) or fax form to 403-253-1744.  
Contact our customer service team at (877) 252-6655.

| Push Chair Specifications | Stroller Only      | Liner Size 1 | Liner Size 2 | Sitter Size 1 | Sitter Size 2 |
|---------------------------|--------------------|--------------|--------------|---------------|---------------|
| Seat to Footwell          | 13"                | 15"          | 15"          | 16"           | 17"           |
| Seat                      | 14"W x 12"D        |              |              |               |               |
| Back Height               | 24"                |              |              |               |               |
| Maximum User Weight       | 90 lbs.            |              |              |               |               |
| Stroller Weight           | 22 lbs.            |              |              |               |               |
| Handle Height to Floor    | 30" - 43"          |              |              |               |               |
| Dimensions when Open      | 45"H x 24"W x 46"D |              |              |               |               |
| Dimensions when Folded    | 16"H x 18"W x 32"D |              |              |               |               |

| Item                           | Part #         | MSRP        |
|--------------------------------|----------------|-------------|
| Special Tomato® EIO Push Chair | _____ 63000000 | \$ 1,455.23 |
| Rain Cover                     | _____ 63030000 | \$ 119.28   |

#### Soft-Touch® Liners Back Cushions

| Size | Cherry  | Lilac  | Aqua  | Dark Gray  | MSRP      |
|------|-------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|-----------|
| LB1  | _____ 72050101                                                                            | _____ 72050102                                                                           | _____ 72050104                                                                           | _____ 72050109                                                                                 | \$ 315.64 |
| LB2  | _____ 72050201                                                                            | _____ 72050202                                                                           | _____ 72050204                                                                           | _____ 72050209                                                                                 | \$ 343.06 |

#### Soft-Touch® Liners Seat Cushions

| Size | Cherry  | Lilac  | Aqua  | Dark Gray  | MSRP      |
|------|--------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|-----------|
| LS1  | _____ 72040101                                                                             | _____ 72040102                                                                            | _____ 72040104                                                                            | _____ 72040109                                                                                  | \$ 315.64 |
| LS2  | _____ 72040201                                                                             | _____ 72040202                                                                            | _____ 72040204                                                                            | _____ 72040209                                                                                  | \$ 343.06 |

#### Soft-Touch® Sitters

| Size   | Lilac  | Aqua  | Dark Gray  | MSRP      |
|--------|-------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|-----------|
| Size 1 | _____ 77010102                                                                            | _____ 77010104                                                                            | _____ 77010109                                                                                  | \$ 660.10 |
| Size 2 | _____ 77010202                                                                            | _____ 77010204                                                                            | _____ 77010209                                                                                  | \$ 868.58 |