

**NON-POWERED**  
Mattress Replacement Systems

**ROHO® SOFFLEX® 2 MATTRESS REPLACEMENT SYSTEM**

The newest therapeutic solution for the superior skin protection you expect from ROHO.

The ROHO SOFFLEX 2 Mattress Replacement System is a superior skin protection support surface for different care settings, costing less than leading powered mattress systems. Utilizing ROHO’s Shape Fitting Technology®, SOFFLEX 2 Mattress System’s interconnected air cell sections may be adjusted separately to match the individual’s shape, thus providing overall skin and soft tissue protection.

**PRODUCT FEATURES**

- ROHO’s Shape Fitting Technology® redistributes forces that can cause peak pressures, shear and friction. The SOFFLEX 2 Mattress is effective for comfort, skin/soft tissue protection and treatment.
- Adjustable – Three sections create a multi-zoned support surface for an individualized fit. The side inflation port allows easy, quick access to adjusting each section.
- Fits most bed frames. Including existing home, hospital or care facility bed frames that are 80” to 84” (203cm to 213cm) in length.
- Cost effective and energy efficient. The SOFFLEX 2 Mattress is non-powered, reusable and costs significantly less than powered mattress systems.
- Lightweight! The SOFFLEX 2 Mattress weighs only 28 lbs. (12.7 kg). Easy to transport and handle during set up and adjustment.
- Convenient Cover. The SOFFLEX2 Mattress’ Cover has a fluid-resistant top and that zips together with the bottom to completely enclose and protect the mattress. The cover is also flame-resistant and the cover top can be wiped clean and removed for disinfecting.

**Size:** SFLX2MRS--STD: 80” x 36” x 8.5” (203cm x 91.5cm x 21.5cm)  
SFLEX2MRS-CUST: 84” x 36” x 8.5” (213.5cm x 91.5cm x 21.5cm)  
SFLX2MRS-3680: 36” x 80” x 8.5” (91.5cm x 203cm x 21.5cm)  
SFLX2MRS-3684: 36” x 84” x 8.5” (91.5cm x 213.5cm x 21.5cm)  
SFLX2MRS-3980: 39” x 80” x 8.5” (99cm x 203cm x 21.5cm)  
SFLX2MRS-3984: 39” x 84” x 8.5” (99cm x 213.5cm x 21.5cm)  
SFLX2MRS-4280: 42” x 80” x 8.5” (106.5cm x 203cm x 21.5cm)  
SFLX2MRS-4284: 42” x 84” x 8.5” (106.5cm x 213.5cm x 21.5cm)  
SFLX2MRS-4880: 48” x 80” x 8.5” (122cm x 203cm x 21.5cm)  
SFLX2MRS-4884: 48” x 84” x 8.5” (122cm x 213.5cm x 21.5cm)

**Approximate Weight:** Approximately 28 lbs. (12.7 kg).

**Weight Limit:** 650 lbs. (295 kg). Mattress should be correctly sized to the bed frame.

**Warranty:** Mattress: 5-Year Limited Warranty, Cover: 18-month Limited Warranty



**U.S. Medicare Code:**  
**E0373 Group 2 Support Surfaces**  
**Available in USA only**



**FOR PEOPLE WHO:**

- are at risk of developing pressure ulcers and/ or tissue injury or who currently have skin/soft tissue breakdown
- have any stage pressure ulcer, skin grafts or recent myocutaneous flap surgery
- need comfort
- need a mattress that provides an individualized fit