

## Non-Padded Transfer Benches

**GUARDIAN®**

- Allows patients to enter and exit tubs safely
- Suction cup feet grip the tub floor for added safety
- Side arm provides extra stability and leverage
- Height-adjustable legs accommodate most tub wall heights



G98309G

### Premium Transfer Bench

Built-in soap dish and hand-held shower slot maximize bathing convenience. Backrest adjusts without tools for right or left facing tubs. Comes in an attractive retail box.



MDS86960R

### Standard Transfer Bench

No tools are required to assemble this durable plastic model that is rustproof. Included with this transfer bench is an optional removable arm that helps with patient balance and a built-in handheld shower holder and soap dish.



G98308A

### Basic Transfer Bench

Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage. Backrest attaches for right or left facing tubs easily without tools.

### Premium Bariatric Transfer Bench

Total of 6 legs and cross brace supporting legs inside tub provide maximum support for larger patients. Direction of the seat back can be reversed easily without tools.





MDS86952XW



MDS86952

### Basic Transfer Bench

Blow molded transfer bench provides ideal support and stability at a more economical price. Direction of the seat back can be reversed easily without tools.

Item Number	Description	Color	Weight Cap Lbs (kg)	Seat Height Inch (cm)	Seat Width Inch (cm)	Seat Depth Inch (cm)	Pkg
<b>Non-Padded Transfer Bench</b>							
G98309G	Premium Non-Padded	White	300 (136)	17-21.5 (43.2-54.6)	31.5 (80)	21 (53.3)	3/cs 
MDS86960R	Standard Composite Plastic	White	350 (159)	17-21 (43.2-53)	17.5 (44.5)	29.5 (75)	2/cs 
G98308A	Non-Padded	Gray	300 (136)	15.5-20.5 (39.4-52.1)	26.5 (67.3)	16.25 (41.3)	2/cs
MDS86952	Basic	White	300 (136)	16.5-21.5 (41.9-54.6)	17 (43.2)	16.5 (41.9)	2/cs
<b>Non-Padded Bariatric Transfer Bench</b>							
MDS86952XW	Premium Bariatric	White	550 (250)	16-20 (40.6-50.8)	17 (43.2)	16 (40.6)	Each

 **Retail Packaging** Ready for Display