



Stingray[™]



Stingray is a superbly stylish buggy with tool free adjustment of postural support and unique swivel action for forward/rearward facing options. The Stingray is available in 2 sizes.

FABRIC



FRAME



180 - 310 mm



180 - 300 mm



430 - 670 mm



max. 40 kg

Foldable

The Stingray is folded quickly and fits into almost any car. Fold the seat and frame units separately or altogether and pack it into your car - soon you will be on your way.



180° turnable seat unit.



Tilt in space and back recline

Perfect seating position due to adjustable seat tilt and back angle adjustment. The seat unit tilts 45° and the back reclines down, allowing your child to rest while out and about.




Safe transportation in vehicles

The Stingray complies with ISO 7176-19, for safe transportation in motor vehicles. The Stingray is approved for transport in vehicles of a user weight of up to 40 kgs.



Accessories



1 Head rest cushion
8740810, 8740820

2 Sculpted head support
872083x

3 Folding top
8740110, 8740120

4 Anatomic head rest
8720862, 8720863

5 Grip handle
872040x

6 Hip supports
87402x0

7 Pommel
872050x

8 Tray
87200x0
8776350

9 Transport fittings
8720902

10 Side supports
872030x

11 Foam back plate
877632x

For added safety and support.

12



13



14



15



16



12 Hip belt
85427-x

13 Oxygen bottle holder
8776295

14 Bunting bag
874013x

15 Vent tray
8776290

16 Shopping basket
8720600