

Drop-Arm Commodes

- Arms swing out of the way for easier lateral transfers on and off the commode
- Transfers are easier for patients using a wheelchair
- Place near the bed with near arm down for patients transferring during the night
- Comes standard with pail, lids and splash guard



G98204

Premium Padded

Thick resilient foam padded seat and back for increased comfort is ideal for long-term use. Hard, chrome-plated steel is rust resistant and offers great durability.



G98202

Standard Non-Padded

Hard, chrome-plated steel frame is rust resistant and offers greater durability.



MDS89668

Basic Economy

Chrome-plated steel frame with easy, single-handed release mechanism for quick transfers. Comes standard with 12-quart pail with cover included.

Item Number	Description	Weight Cap Lbs (kg)	Seat Height Range Inch (cm)	Seat Depth Inch (cm)	Width Between Arms Inch (cm)	Overall Width Inch (cm)	Pkg
G98204	Premium Padded	350 (159)	20-25 (50.8-63.5)	18-20 (45.7-50.8)	18.5 (47)	22 (55.9)	Each
G98202	Standard Non-Padded	300 (136)	18.25-23.25 (46.4-59.1)	18-20 (45.7-50.8)	18.5 (47)	22 (55.9)	Each
MDS89668	Basic Economy	250 (113)	17-21 (43.2-53.3)	15 (38.1)	17.5 (44.5)	25.5 (64.8)	Each