

TILT REHAB SHOWER COMMODES

Traditionally, rehab shower commodes have been constructed of steel or stainless steel. These two materials had their drawbacks. Steel rusts, stainless steel is expensive and both are heavy. Ezee Life™ commodes are constructed of high quality aluminum that won't rust, is lighter than steel and lower in cost than stainless steel.





Model # | 155 (15" seat)
Model # | 190 (18" seat)
Model # | 195 (20" seat)



Model # | 155-22 (15" seat, 22" wheels)
Model # | 190-24 (18" seat, 24" wheels)
Model # | 195-24 (20" seat, 24" wheels)



Optional 4 Way Seat



Anti tippers
on all models



Height adjustable
on all models



Open Front Seat & Insert



Standard on all sizes

Closed Front Seat & Top



Available on 190 & 190-24

4-Way Open Seat



Available on 190 & 190-24

Removable snap on seat (all models)



STANDARD FEATURES:

- Lightweight aluminum **NO RUST**
- Footrest, adjustable height & removable
- Removable insert cushion for showering or toileting
- Quick Release axles (self propelled models)
- Adjustable headrest up/down/forward/back
- Four 5" locking casters
- Slide on pail, removable
- Flip back arms
- Maintenance free tires
- Calf strap included
- Molded handrims
- Super soft seat
- Easier rolling
- Seat belt
- All adjustments are tool-less
- 2 year Warranty

SPECIFICATIONS	155 Tilt	190 Tilt	195 Wide Tilt
Frame Material	Aluminum	Aluminum	Aluminum
Weight Capacity	286 lbs (130 kgs)	300 lbs (136 kgs)	325 lbs (148 kgs)
Rear Wheels	5" locking casters	5" locking casters	5" locking casters
Optional Rear Wheels	22" (155-22)	24" (190-24)	24" (195-24)
Seat Width	15"	18"	20"
Seat Depth	15"	17"	17"
Foot Rest Hanger Angle	60°	60°	60°
Foot Rest Length	14" to 20"	14" to 20"	14" to 20"
Back Height	18"	18"	18"
Product Weight	32.5 lbs (14.8 kgs)	36 lbs (16.4 kgs)	38 lbs (17.3 kgs)
Arm Height	11" (from seat)	11" (from seat)	11" (from seat)
Arm Length	13"	13"	13"
Arm Pad Length	9"	9"	9"
Tilt Range Posterior	20°	20°	20°
Rear Internal Width	15"	15"	17"
Seat Heights (Standard)	21", 22" or 23"	21", 22" or 23"	21", 22" or 23"
Seat Heights (Tall)	22.9", 23.9", 24.9"	22.9", 23.9", 24.9"	22.9", 23.9", 24.9"
Seat Clearance Height ²	19", 20", 21"	19", 20", 21"	19", 20", 21"
Seat Clearance Height ²	21", 22", 23"	21", 22", 23"	21", 22", 23"
DIMENSIONS WITH 5" REAR CASTERS	155 Tilt	190 Tilt	195 Wide Tilt
Overall Width	19"	22"	24"
Overall Length ¹	40"	40"	40"
DIMENSIONS WITH 24" REAR WHEELS	155 Tilt	190 Tilt	195 Wide Tilt
Overall Width	23"	26"	28"
Overall Length ¹	45"	45"	45"

- 1) Length with all parts extended (wheels highest height, anti-tip bars fully extended, footrests in place)
 2) Seat clearance is the highest toilet bowl height the commode can be positioned over.