

Clever-Lite LS, Adult Walker

HCPCS: E0143



FIGURE A



FIGURE B

- Rear glide push-down brakes (Figure A)
- Sturdy 1" diameter anodized, extruded aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame
- Easily folds with patented dual-lever side-paddle release (Figure B)
- Soft, flexible backrest provides comfort and security while seated
- 5" fixed front wheels

PRODUCTS	Description
10233	Clever-Lite LS Rollator, w/ 5",1/cs
SPECIFICATIONS	Description
Carton Shipping Weight	15 lbs.
Handle Height	29.5" - 39"
Seat Dimensions	18.75" (W) x 11" (D)
Seat Height	19" - 25.5"
Warranty	Limited Lifetime
Weight	11 lbs.
Weight Capacity	300 lbs
Width	18.5" inside hand grip 21" inside back legs 27" Overall
PARTS	Description
10230SN	Blow Molded Seat,10230,1ea, EA
10109	Walker Wheels,Univrsl 5",1pr/bx, RTL, BX
10230SF	Seat Frame,10230,1ea, EA
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
10233RL	Rear Legs With Glide Brake,10233,1pair, PR
10325	WalkerTipGrey,10233,1/ea, EA
10325B	WalkerTipBlack,11/8". 1/ea, EA
12002PB	Push Button for Bath Bench,Walker,1/ea, EA
950300169901	Anti Rattler, 795/10215/10216, 2/pk, EA