

Folding Universal Sliding Transfer Bench



FIGURE A



FIGURE B



FIGURE C



RETAIL BOX

- Transfer Bench slides left and right along frame for easy entry into and out of tub
- Seat height adjusts in ½" increments
- Folds flat for transport and storage (Figure A)
- Comes standard with removable soap dish (Figure B)
- Comes with safety net to catch dropped personal items and prevent bending
- Comes with seat belt for added safety
- Retail packaged (Figure D)

PRODUCTS

Description

RTL12075	Sliding Shower Chair w/Soap Box, 2/cs
----------	---------------------------------------

SPECIFICATIONS

Description

Seat Dimensions	22.5" (W) x 18.5" (D)
Height	21" - 25"
Weight	16 lbs.
Carton Shipping Weight	38 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12075-A	Replacement Handle, 1/ea
12075B	Replacement Back, 1/ea
12075S-SD	Replacement Soap Dish, 1/ea