

Plastic Tub Rail

HCPCS: E0246



FIGURE A



FIGURE B

- Attractive, ergonomically-designed tub rail provides a comfortable, safe grip
- Handle is perpendicular to tub
- 16.5" in height
- Tool-free, easy to install
- Width adjusts from 3" - 6.25"
- Retail packaged
- Made using recycled materials

PRODUCTS	Description
12044	Tub Rail,16 1/2"Plastic,1/bx, RTL

SPECIFICATIONS	Description
Adjustable Clamp	3" (smallest) 6.25" (widest)
Carton Shipping Weight	7 lbs.
Dimensions	15" (W) x 12.5" (H) x 4" (D)
Height	16.5"
Warranty	Limited Lifetime
Weight	5 lbs.
Weight Capacity	300 lbs