

Premium Plastic, Raised, Elongated Toilet Seat with Lock



- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms are 21" apart and can be removed or added as needed
- Wide opening in front and back for personal hygiene
- Heavy-duty molded plastic construction provides additional strength and durability
- Standard locking mechanism ensures safety
- Lightweight and portable
- No tools required for installation
- Easy-to-clean
- Fits most elongated toilets

PRODUCTS

Description

12013	With Arms, 1/cs
-------	-----------------

SPECIFICATIONS

Description

Seat Dimensions	15.5" (W) x 17" (D) x 5" (H)
Width Between Arms	21"
Weight	4.5 lbs.
Carton Shipping Weight	6.5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime