

BANANA BOARD PRIME

The Banana Board Prime is a comfortable, lightweight transfer board, designed to aid seated transfers of individuals between two surfaces, such as from bed to wheelchair, wheelchair to chair and wheelchair to car.

Features:

- Can be used for both assisted and non-assisted transfers
- Supports independent movement
- Curved board shape allows side-to-side and straight transfers
- Tapered edges to ensure easy insertion and removal
- Non-slip pads to ensure safe movement
- Smooth moulded surface
- Ultra-strong and long lasting Ideal for most users
- Ideal for both adults and paediatric users
- Lightweight and easy to carry
- Easy to clean
- Maximum weight limit 200kg



Specifications:

Measurements: 62cm long x 26 cm wide

Thickness: 1cm

Weight: 0.8kg

Weight Limit: 200kg