



**BRODA<sup>®</sup>**  
*Comfort is Essential*

**TRANQUILLE**  
AUTO-LOCKING GLIDER



## NURTURING CONTENTMENT AND INDEPENDENCE

A stationary chair, the Auto Locking Glider provides gentle gliding for use in patient rooms, common areas and memory care.



- 1 HEIGHT ADJUSTABLE HEADREST (STANDARD ON HIGH-BACK)
- 2 HEIGHT ADJUSTABLE ARMRESTS
- 3 AUTO-LOCKING MECHANISM
- 4 CAREGIVER LOCK (MAKES CHAIR STATIONARY)
- 5 REAR TRANSPORT WHEELS

### PATIENT BENEFITS



Gentle gliding motion helps calm and soothe agitation.



Auto-locking features aids in patient safety and fall prevention.



Gliding motion provides stimulation for active patients helping reduce the need to wander.



Rear transport wheels allow the glider to be moved from room to room with ease.

### CAREGIVER AND FACILITY BENEFITS

- Smooth gliding motion provides relaxation
- Gentle movement helps prevent wandering
- Auto-locking features aids in patient safety and fall prevention
- Provides stimulation and contentment - Excellent for memory care
- BRODA's Comfort tension seating provides seating comfort and pressure redistribution



## COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue



CTS Strapping flexes for individualized fit & superior comfort.

## SPECIFICATIONS


Frame	16 Gauge tubular steel frame
Seat & Backrest	BRODA Comfort Tension Seating®
Wheels	Rear transport wheel, wheels for non-patient transport
Weight Capacity	250 lbs. (114 kg)
Standard Padding	Seat and back pad vinyl cover
HCPCS	A9270

Model	Seat width	Seat depth	Seat height	Back height
Small low-back 100-10AL	18"	17"	16"	26"
Low-back 100-15AL	20"	17"	18"	26"
High-back 100-20AL	20"	17"	18"	32"


## COMMON ACCESSORIES

- Padded armrests
- Depth & height adjustable armrest
- Full or half tray
- Clear tray
- Contoured neck rest (high back glider only)

## UPHOLSTERY OPTIONS



RENAISSANCE  
BURGUNDY



RENAISSANCE  
CHARCOAL



RENAISSANCE  
METEOR



SILVERTEX  
SAPPHIRE BLUE

## FRAME OPTIONS



GREY  
HAMMER



CHESTNUT  
BROWN





 **MITY**, INC  
HEALTHCARE DIVISION  
[info@brodaseating.com](mailto:info@brodaseating.com)  
p. 1-844-552-7632

[Brodaseating.com](http://Brodaseating.com)  
©Broda USA, INC.  
8/16