



STANDARD MODEL FEATURES

- Support users up to 66" in height and up to 130 lb in weight
- Low profile sides for easy transfer of users
- Durable, machine-washable padded cover
- Memory foam back/seat padding reduces sores and fatigue
- Soft-padded headwings
- Latex and silicone free
- Crash tested to FMVSS #213

PLUS MODELS FEATURE

1. Swing-Away Padded Trunk Supports
2. Swing-Away Padded Hip Supports
3. Swing-Away Abductor



UPHOLSTERY



Convertible Pink
(CP)



Dune Buggy Beige
(DB)



Go Kart Green
(GG)



Roadster Red
(RR)



Speedway Gray
(SG)

ACCESSORIES



Incontinence Covers



3-inch Extensor
Thrust Wedge



Buckle Guard



Retainer Clip Guard

ORDER FORM FAX TO 905-787-8550 OR SUBMIT ONLINE

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

BASE UNIT	CSS-2400S	CSS-2400M	CSS-2400L	CSS-2401M	CSS-2401L
Product Size	Standard w/o Seat Depth Extender	Plus-M w/o Seat Depth Extender	Plus-L w/o Seat Depth Extender	Plus-M <u>with</u> Seat Depth Extender	Plus-L <u>with</u> Seat Depth Extender
Weight Capacity	25 - 130 lb	25 - 130 lb	25 - 130 lb	25 - 130 lb	25 - 130 lb
Max. Height Capacity	66"	66"	66"	66"	66"
Harness Height	11", 13.5", 16", 18.5", 20.5", 23.5" (min. 10")				
Shoulder Width	18"	18"	18"	18"	18"
Back Height	28"	28"	28"	28"	28"
Seat Depth	12"	12"	12"	16"	16"
Inside Seat Width at Front	17"	17"	17"	19"	19"
Inside Seat Width at Rear	15"	15"	15"	15"	15"
Overall Width	21.5"	21.5"	21.5"	21.5"	21.5"
Seat-to-Back Angle	98°	98°	98°	98°	98°
Product Weight	27 lb	29 lb	29 lb	33 lb	33 lb

PLUS PACKAGE INCLUDES:

Select to have PLUS components itemized on quote.
(quote will show standard car seat item number (CSS-2400S) + components below)

Seat Depth Extender	N/A	N/A	N/A	☒ CSS-R2425	☒ CSS-R2425
Abductor	N/A	☒ CSS-R2415	☒ CSS-R2415	☒ CSS-R2419	☒ CSS-R2419
Hip Bracket	N/A	☒ CSS-R2473H Width: 10" - 14"	☒ CSS-R2471H Width: 13" - 16"	☒ CSS-R2473H Width: 10" - 14"	☒ CSS-R2471H Width: 13" - 16"
Hip Pad	N/A	☒ CSS-R2481M Pad: 4" x 8"	☒ CSS-R2481L Pad: 4" x 12"	☒ CSS-R2481M Pad: 4" x 8"	☒ CSS-R2481L Pad: 4" x 12"
Trunk Bracket	N/A	☒ CSS-R2474T Width: 6" - 12"	☒ CSS-R2473T Width: 10" - 14"	☒ CSS-R2474T Width: 6" - 12"	☒ CSS-R2473T Width: 10" - 14"
Trunk Pad	N/A	☒ CSS-R2482M Pad: 3" x 5"	☒ CSS-R2482L Pad: 4" x 6"	☒ CSS-R2482M Pad: 3" x 5"	☒ CSS-R2482L Pad: 4" x 6"

UPHOLSTERY

CP 	DB 	GG 	RR 	SG 
Convertible Pink with Gray Mesh	Dune Buggy Beige with Gray Mesh	Go Kart Green with Gray Mesh	Roadster Red with Black Mesh	Speedway Gray with Gray Mesh

ACCESSORIES

Incontinence Cover Kit	CSS-A2405S	CSS-A2405M	CSS-A2405L	CSS-A2405M	CSS-A2405L
Seat Depth Extender	CSS-R2425	(increases 4" to seat depth, comes standard on CSS-2401 packages)			
Extensor Thrust Wedge	CSS-A2455	(raises 3" on the front)			
Tether Strap Kit	CSS-A2030	(one kit included with car seat)			
Bolt-in Hardware Kit	CSS-A2031	(for vehicle manufactured before 2002)			
Retainer Clip Guard	CSZ-A2085				
Buckle Guard	ZZ-A2075				