

R82 Flamingo | Flamingo High-Low

Order Form



Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:


Flamingo Standard Features:


- Size 1, 2, 3, 4
- 4 White Locking Castors
- Aluminum Frame
- 100% Waterproof Upholstery for Seat and Back
- Upholstery Available in Blue
- Height Adjustable and Foldable Frame
- Tilt in Space (+5° to -35°)
- Reclining Back
- Active Back (No Sides)
- Splash Guard - High Model



Flamingo High-Low Standard Features:

- Size 1, 2, 3, 4
- 4 White Locking Castors
- Foot Support
- Stainless Steel Frame
- 100% Waterproof Upholstery for Seat and Back
- Upholstery Available in Blue
- Height Adjustable Frame
- Tilt in Space
- Reclining Back
- Active Back (No Sides)
- Splash Guard - High Model

<u>Flamingo Seat and Back</u>		<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A) 	<input type="checkbox"/> Flamingo Seat and Back, Blue (A)	1		880101-61
	<input type="checkbox"/> Flamingo Seat and Back, Blue (A)	2		880102-61
	<input type="checkbox"/> Flamingo Seat and Back, Blue (A)	3		880103-61
	<input type="checkbox"/> Flamingo Seat and Back, Blue (A)	4		880104-61



<u>Flamingo Frame</u>		<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A) 	<input type="checkbox"/> Flamingo Frame w/out Footrest (A)	1		880201-NL
	<input type="checkbox"/> Flamingo Frame w/out Footrest (A)	2		880202-NL
	<input type="checkbox"/> Flamingo Frame w/out Footrest (A)	3		880203-NL
	<input type="checkbox"/> Flamingo Frame w/out Footrest (A)	4		880204-NL







<u>Flamingo High-Low Bath Frame</u>		<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A) 	<input type="checkbox"/> High-Low Bath Frame w/ Footrest (A)	1		880411
	<input type="checkbox"/> High-Low Bath Frame w/ Footrest (A)	2		880412
	<input type="checkbox"/> High-Low Bath Frame w/ Footrest (A)	3		880413
	<input type="checkbox"/> High-Low Bath Frame w/ Footrest (A)	4		880414





<u>Back Option</u>		<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)  (B) 	<input type="checkbox"/> Back Extension (H: 4" x L: 9") (A)*	1, 2		880118-61
	<input type="checkbox"/> Back Extension (H: 4" x L: 9" - 10.2") (A)*	3, 4		880119-61
	<input type="checkbox"/> Mounting Rod for Back Extension	7.9"		880120
	<input type="checkbox"/> Back Reinforcement Kit (B)	3, 4		880155



(A)* Requires Mounting Rod for Back Extension






		<u>Pelvic Positioning</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)		(B) 	<input type="checkbox"/> Upholstery, Seat w/ Wide Sides, Blue (A)*	2	880166-61
			<input type="checkbox"/> Upholstery, Seat w/ Wide Sides, Blue (A)*	3	880167-61
			<input type="checkbox"/> Upholstery, Seat w/ Wide Sides, Blue (A)*	4	880168-61
			<input type="checkbox"/> Seat Infill for Seat, Blue (B)	1	880185-61
			<input type="checkbox"/> Seat Infill for Seat, Blue (B)	2	880186-61
			<input type="checkbox"/> Seat Infill for Seat, Blue (B)	3	880187-61
			<input type="checkbox"/> Seat Infill for Seat, Blue (B)	4	880188-61
			<input type="checkbox"/> Seat Belt w/ Plastic Lock	One Size	85100
(A)* Does not include Splashguard. Reduces sz 2 & 3 width by 2", and sz 3 by 3.5"					



		<u>Foot & Leg Positioning</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)		(B) 	<input type="checkbox"/> Flamingo Foot Support (A)	1	880211
			<input type="checkbox"/> Flamingo Foot Support (A)	2	880212
			<input type="checkbox"/> Flamingo Foot Support (A)	3	880213
			<input type="checkbox"/> Flamingo Foot Support (A)	4	880214
			<input type="checkbox"/> Calf Support (H: 3.9" x L: 41.7") (B)	One Size	880199
			<input type="checkbox"/> Footrest Straps, Pair		85096

		<u>Trunk Positioning</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)		(B) 	<input type="checkbox"/> Bar, Straight, Stainless (A)**	3¾"	82569
			<input type="checkbox"/> Bar, Straight, Stainless (A)**	6¾"	82648
			<input type="checkbox"/> Cushion, Curved, Each (B)**	A: 3¾", B: 7½"	6175407
			<input type="checkbox"/> Cushion, Curved, Each (B)**	A: 2½", B: 6½"	6175408
(C)		(D) 	<input type="checkbox"/> Cushion, Curved, Each (B)**	A: 2½", B: 5¾"	6175409
			<input type="checkbox"/> Protection Pads, Each		89151
			<input type="checkbox"/> Swing-Away Bracket, Right, Each (C)**	A: 1¾", B: 4¾"	83300-R
			<input type="checkbox"/> Swing-Away Bracket, Left, Each (C)**	A: 1¾", B: 4¾"	83300-L
			<input type="checkbox"/> Swing-Away Bracket, Right, Offset, Each (D)	A: 2¼", B: 5"	83500-R
			<input type="checkbox"/> Swing-Away Bracket, Left, Offset, Each (D)	A: 2¼", B: 5"	83500-L
(E)		(F) 	<input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)*	1	880111-61
			<input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)*	2	880112-61
			<input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)*	3	880113-61
			<input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)*	4	880114-61
			<input type="checkbox"/> Vest, Water-Proof Mesh (F)	0	89111-040
			<input type="checkbox"/> Vest, Water-Proof Mesh (F)	¾"	89111-07540
			<input type="checkbox"/> Vest, Water-Proof Mesh (F)	1	89111-140
			<input type="checkbox"/> Vest, Water-Proof Mesh (F)	2	89111-240
			<input type="checkbox"/> H-Harness w/ Plastic Lock		85105
(A), (B), (C)** Required to complete Swing-Away Lateral Assy					
(E)* Cannot be combined w/ Swing-Away Lateral Supports					

		<u>Head Support</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)		(B) 	<input type="checkbox"/> Head Support No. 1, Blue PU (A)		6071000
			<input type="checkbox"/> Head Support No. 2, Blue PU (B)		6072000
(C)		(D) 	<input type="checkbox"/> Head Support No. 3, Blue PU (C)		6073000
			<input type="checkbox"/> Head Support No. 4, Blue PU (D)		6074000
			<input type="checkbox"/> Turtle Bar System, Short (E)	5.1"	98619-25
			<input type="checkbox"/> Turtle Bar System, Long (E)	9.8"	98619-26

		<u>Frame Options</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)			<input type="checkbox"/> Commode Pan w/ Holder for Flamingo (A)	1	880121
			<input type="checkbox"/> Commode Pan w/ Holder for Flamingo (A)	2	880122
			<input type="checkbox"/> Commode Pan w/ Holder for Flamingo High-low (A)	1	880121-HL
			<input type="checkbox"/> Commode Pan w/ Holder for Flamingo High-low (A)	2	880122-HL
			<input type="checkbox"/> Commode Pan w/ Holder for Flamingo & High-Low (A)	3	880123
			<input type="checkbox"/> Commode Pan w/ Holder for Flamingo & High-Low (A)	4	880124

Frame Options Continued		Size	HCPCS	Part No.
(B) 	<input type="checkbox"/> Potty Ring w/ Holders (B)*	1		880121-PR
	<input type="checkbox"/> Potty Ring w/ Holders (B)*	2		880122-PR
	<input type="checkbox"/> Potty Ring w/ Holders (B)*	3		880123-PR
	<input type="checkbox"/> Potty Ring w/ Holders (B)*	4		880124-PR
(C) 	<input type="checkbox"/> Push Brace (C)	One Size		880110
(D) 	<input type="checkbox"/> Push Handle for Flamingo High-Low (D)	1, 2		880110-HL
	<input type="checkbox"/> Push Handle for Flamingo High-Low (D)	3, 4		880110-HLL
(E) 	<input type="checkbox"/> Single Post Armrest (L: 9"), Set (E)	One Size		880152
	<input type="checkbox"/> Suction Cups for Seat (F)	One Size		880150
(F) 				
(B)* Potty ring required for use over the commode				
(E)* Armrest cannot be used together w/ Sides for Back - 88011X-X				

Accessories		Size	HCPCS	Part No.
(A) 	<input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A)	1		880131-61
	<input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A)	2		880132-61
	<input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A)	3		880133-61
	<input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A)	4		880134-61
(B) 	<input type="checkbox"/> Splashguard, Low Model, Blue (B)	1		880195-61
	<input type="checkbox"/> Splashguard, Low Model, Blue (B)	2, 3, 4		880196-61

IMPORTANT: No returns or exchanges of hygiene products

For shipping rates, please call 1-844-876-6245

PO Box 1739 Matthews, NC 28106 • 1-844-876-6245 • Fax: 704-882-0751


Flamingo Configuration

- ☐ Flamingo/Flamingo High-low
- ☐ Head support
- ☐ Push handle/push brace
- ☐ Side supports
- ☐ Arm rest/tray
- ☐ Commode pan/potty ring
- ☐ Vests and belts



A (Back height) _____
 B (Arm height) _____
 C (Lower leg length) _____
 D (Seat depth) _____
 E (Chest width) _____
 F (Seat width) _____
 G (Height from seat to top of head) _____
 H (Shoulder width) _____
 Age _____
 Weight _____

Measurements



	Size 1 Inches	Size 2 Inches	Size 3 Inches	Size 4 Inches
Seat width (A)	9¼"	11"	12½"	13¾"
Seat depth (B)	9¾"	11¼"	13¾"	15½"
Back height (C)	13¾"	14¾"	16½"	18¼"
Lower leg length (D)	18" - 26"	25" - 35"	30" - 40"	30" - 40"
-	7" - 10"	9¾" - 13½"	11¾" - 15½"	11¾" - 15½"
Seat height (E)	50" - 60"	50" - 60"	50" - 60"	50" - 60"
-	19½" - 23½"	19½" - 23½"	19½" - 23½"	19½" - 23½"
Back width (F)	7¾"	8½"	10¼"	11½"
Width, frame (G)	17"	18"	19¼"	21"
Length, frame (H)	26"	26"	26"	26"
Height, frame w/seat (I)	33½"	35"	36½"	37"
Seat height, folded	9½"	10"	10½"	10½"
Seat length, folded	14½"	15½"	17"	18½"
Frame height, folded	8"	8"	8"	8"
Frame length, folded	29"	29"	29P	29"
Weight	24 lbs.	29 lbs.	31 lbs.	33 lbs.
Max. Load/user weight	77 lbs.	110 lbs.	132 lbs.	154 lbs.
Angling				
Seat		-15°-+35°	Back	-90°-30°
Foot support bar			90°-60°	

Flamingo High-Low Configuration

- ☐ Flamingo/Flamingo High-low
- ☐ Head support
- ☐ Push handle/push brace
- ☐ Side supports
- ☐ Arm rest/tray
- ☐ Commode pan/potty ring
- ☐ Vests and belts



A (Back height) _____

B (Arm height) _____

C (Lower leg length) _____

D (Seat depth) _____

E (Chest width) _____

F (Seat width) _____


G (Height from seat to top of head) _____

H (Shoulder width) _____

Age _____

Weight _____

Measurements



	Size 1 Inches	Size 2 Inches	Size 3 Inches	Size 4 Inches
Seat width (A)	9¼"	11"	12½"	13¾"
Seat depth (B)	9¾"	11¼"	13¾"	15½"
Back height (C)	13¾"	14¾"	16½"	18¾"
Lower leg length (D)	7¼" - 9¾"	9" - 13¼"	11¾" - 16½"	11¾" - 16½"
Seat height above floor (E)	20¾" - 28¾"	20¾" - 28¾"	20¾" - 28¾"	20¾" - 28¾"
Back width (F)	7¼"	8½"	10¼"	11"
Frame width (G)	19½"	19½"	22¼"	22¼"
Frame length (H)	25"	25"	29½"	29½"
Total height (I)	33½" - 42"	35¾" - 43¼"	36¾" - 45¾"	39¾" - 46¾"
Seat, total width	13¾"	14¾"	16¾"	18¾"
Seat, total depth	10¼"	12¼"	14"	16½"
Seat, total height	15½"	16¾"	18¾"	19¾"
Back angling	-90°-+30°	-90°-+30°	-90°-+30°	-90°-+30°
Tilt in space	-5°-+20° *	-5°-+20°	-5°-+20°	-5°-+20°
Foot support angle	90°-75°	90°-75°	90°-75°	90°-75°
Weight, total	31 lbs.	32 lbs.	36½ lbs.	37½ lbs.
Weight, frame	25½ lbs.	25½ lbs.	27½ lbs.	27½ lbs.
Weight, seat	5½ lbs.	6½ lbs.	9 lbs.	10 lbs.
Max. load, total	77 lbs.	110 lbs.	132 lbs.	154 lbs.
Max. load, seat	77 lbs.	110 lbs.	132 lbs.	154 lbs.
Max. load, operating force	77 lbs.	77 lbs.	121 lbs.	121 lbs.

* The -5° tilt is obtainable at 55 cm (19½") seat height above floor