


# Wave dimensions

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 <p><b>Key user dimension: height</b> Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth.</p> <p><b>Important:</b> Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.</p>			
Item dimensions (inches)	Z210 small	Z220 medium	Z230 large
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½ – 18½	12½ – 18½	12½ – 18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Base width	18	18	21
Base length (legs down) with legs flat	24 29½	24 29½	24 29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max. tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lb)	75	150	250
Back and seat angle adjustments	0°, 15°, 30°, 45°, 60°, 75°		



All three sizes are available in blue and pink.

# 2021 order form

Use dimension chart to select appropriate size.

## Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

### Fabric without calf rest

Blue Z268  
Pink Z269

(Includes seat, back, cover and chest strap)



### Fabric with calf rest

Blue Z266  
Pink Z267

(Includes seat, back, cover and chest strap)



		Z210 small	Z220 medium	Z230 large
* Chair		Z211 <input type="checkbox"/>	Z221 <input type="checkbox"/>	Z231 <input type="checkbox"/>
* Fabric cover	CHOOSE			
	Blue	Z244 <input type="checkbox"/>	Z246 <input type="checkbox"/>	Z248 <input type="checkbox"/>
	Pink	Z245 <input type="checkbox"/>	Z247 <input type="checkbox"/>	Z249 <input type="checkbox"/>
* Chest strap	CHOOSE			
	With lateral positioning	Z212 <input type="checkbox"/>	Z222 <input type="checkbox"/>	Z232 <input type="checkbox"/>
	Without lateral positioning	Z213 <input type="checkbox"/>	Z223 <input type="checkbox"/>	Z233 <input type="checkbox"/>
Head blocks		Z214 <input type="checkbox"/>	Z214 <input type="checkbox"/>	Z234 <input type="checkbox"/>
Leg straps		Z215 <input type="checkbox"/>	Z225 <input type="checkbox"/>	Z235 <input type="checkbox"/>
Calf rest (Calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 <input type="checkbox"/>	Z226 <input type="checkbox"/>	Z236 <input type="checkbox"/>
Tub stand		Z217 <input type="checkbox"/>	Z217 <input type="checkbox"/>	Z217 <input type="checkbox"/>
Shower stand		Z218 <input type="checkbox"/>	Z218 <input type="checkbox"/>	Z218 <input type="checkbox"/>
Tub transfer base (calf rest recommended)	CHOOSE			
	Counterclockwise into tub	Z227 <input type="checkbox"/>	Z227 <input type="checkbox"/>	Z227 <input type="checkbox"/>
	Clockwise into tub	Z228 <input type="checkbox"/>	Z228 <input type="checkbox"/>	Z228 <input type="checkbox"/>
Legs for tub transfer base	CHOOSE			
	Standard	Z224 <input type="checkbox"/>	Z224 <input type="checkbox"/>	Z224 <input type="checkbox"/>
	Long (tub edge over 20" high)	Z238 <input type="checkbox"/>	Z238 <input type="checkbox"/>	Z238 <input type="checkbox"/>
Additional head blocks		Z214 <input type="checkbox"/>	Z214 <input type="checkbox"/>	Z234 <input type="checkbox"/>
Additional fabric cover		Z219 <input type="checkbox"/>	Z229 <input type="checkbox"/>	Z239 <input type="checkbox"/>
Additional chest strap	CHOOSE			
	With lateral positioning	Z212 <input type="checkbox"/>	Z222 <input type="checkbox"/>	Z232 <input type="checkbox"/>
	Without lateral positioning	Z213 <input type="checkbox"/>	Z223 <input type="checkbox"/>	Z233 <input type="checkbox"/>
		<b>Z210</b>	<b>Z220</b>	<b>Z230</b>

\* Required