

**Rebel Wheelchair**  
Single Axle



A



B

Item #	Seat Width	Armrest	Product Front Riggings	UOM
RTLREB18DDA-SF	18"	Removable Desk Arms	Swing-away Footrest	1/cs

- Quick-release rear wheels allow for easy storage or transport. (Figure A)
- Fold-down back and folding frame allows transport in almost any trunk. (Figure B)
- Easily removable swing-away footrests are also convenient for storage.
- Tool-free removable front casters.
- Desk-style armrests make transfers easy and accessibility better.
- Lightweight and extremely maneuverable, easy to navigate tight areas.
- Large accommodating seating surface.
- Carton size makes the Rebel perfect for retail.

#### specifications

Seat Width	18"	Overall Length (with riggings)	43"
Overall Width	25.5" (open)   10.5" (closed)	Weight	37 lbs.
Overall Height	35"	Carton Shipping Weight	42 lbs.
Seat-To-Floor Height	19.5"	Weight Limit	300 lbs.
Back Height	17"	Warranty	Limited Lifetime
Seat Dimensions	18" (W) x 16" (D)		

#### accessories

##### Seat Belts

STDS850	Seat Belts, Auto-Clasp Type Closure	
STDS851	Seat Belts, Hook and Loop Fastener Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric	1/ea

##### Universal Oxygen "E" Cylinder Holder

STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea

##### Wheel Lock Extension

STDS801GT	Wheel Lock Extension	1/ea
-----------	----------------------	------

##### Anti-Fold Bar

STDS806	Anti-Fold Bar	1/ea
---------	---------------	------

##### Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

##### Heel Loops

STDS831	Heel Loops	
---------	------------	--

##### Cup Holder

STDS1040S	Universal Cup Holder	1/ea
-----------	----------------------	------