















Half-Length Bed Rail

Tool-Free Adjustable Width, with Brown-Vein Finish

Item #	Color	UOM	
15208BV	Brown-Vein	1 pr/cs	

- Half rail was designed to comply as part of the FDA Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment.
- Tool-free width adjustment allows the rail to be adjusted in 2" increments from 36" – 42". (Figure A - B)
- Ideal for protection when bed is elevated.
- Chrome finish on second and fifth bar allows for scratch-proof movement of adjustment bracket.
- Fits most beds with steel frame. (Figure C)
- Attractive brown-vein finish is easy to maintain and less institutional looking. (Figure D)
- Constructed of 1" steel.
- Not for use with "Pan-Style" bed surfaces.
- Easy to install.

specifications	
Length Adjustment	36", 38", 40", 42"
Warranty	Limited Lifetime

Foot Stool

Item #	Color	UOM
13030-1SV	Silver Vein	1/cs
13030-2SV	Silver Vein	2/cs

- Steel-welded silver vein construction. (Figure A)
- Non-skid ribbed rubber platform.
- Reinforced rubber tips.

specifications	
Step Surface	14" (W) x 11" (D)
Step Height	9"
Overall Dimensions	17" (W) x 13" (D)
Carton Shipping Weight (1 2)	7 lbs. 14 lbs.
Weight	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime