

Airgo ProCare IC (Infection Control) Wheelchair

HCPCS:



FIGURE A



FIGURE B



- Designed based on feedback from healthcare professionals, our Airgo ProCare IC wheelchair exceeds institutional standards of quality and features touch points which are easy to both clean and disinfect
- It's streamline design along with mag wheels combine to make it lighter and more durable
- Features dual axle positions that allow for varying seat-to-floor heights
- Heavy-duty chrome-plated steel frame
- Solid tires with adjustable brakes and sealed bearing assemblies
- Ergonomic handles, and flame-retardant black naugahyde upholstery
- Institutional tested durability & proven quality
- Proprietary contact points are easy to clean & disinfect
- Latex free

PRODUCTS	Description
700-622	Airgo ProCare IC (Infection Control) Wheelchair, EA
700-625	Airgo ProCare IC (Infection Control) Wheelchair, EA
700-633	Airgo ProCare IC (Infection Control) Wheelchair, EA

SPECIFICATIONS	Description
Folded Width	12" to 13"
Handle Height	36"
Seat Depth	16" to 18"
Seat Height	20" or 18"
Seat Width	16" 18" 20" 22"
Warranty	Frame and Cross-brace 3 Years, All Other Metal Parts 6 Months, Upholstery, Plastic and Rubber 3 Months
Weight	49 lbs. to 66 lbs. net
Weight Capacity	300 lbs