



These foams are lightweight, yet strong and flexible

Closed cell, water-repellant and mould and mildew resistant

They can be laminated, die-cut, printed and the thermo foams may be dry heated and formed

Excellent for splints, orthotics and prosthetics

2.5lb Density

L-200 Blue - Half Sheet

#7638 – ¼" x 24" x 48"

#7630 – ½" x 24" x 48"

#7631 - 1" x 24" x 48"

#7632 - 2" x 24" x 48"

#7633 - 3" x 24" x 48"

L-200 Blue - Full Sheet

#7644 – ¼" x 48" x 72" (will be cut 2 pieces of 48" x 36")

#7640 – ½" x 48" x 72" (will be cut 2 pieces of 48" x 36")

#7641 - 1" x 48" x 72" (will be cut 2 pieces of 48" x 36")

#7642 - 2" x 48" x 72" (will be cut 2 pieces of 48" x 36")

#7643 - 3" x 48" x 72" (will be cut 2 pieces of 48" x 36")