

Set n' Go Height Adjustable Rollator

The Lumex Set n' Go Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one. Designed with multiple height settings on the frame, the Set n' Go allows the seat height to be adjusted from a hemi 18" to a standard 22" height in 1" increments. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated.

- Easy height adjustment on the frame allows seat height to be adjusted from 18" to 22" in 1" increments
- Comfortable, wide, seamless padded seat with zippered pouch under seat for security of personal items
- Also includes black pouch for storage of larger items
- Ergonomic handgrips offer a comfortable grip
- Adjustable handle height with knobs to secure in place
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Removable padded backbar with easy release pull knobs
- 6" wheel for both indoor and outdoor use
- 300 lb maximum weight capacity
- Designed to meet the requirements of HCPCS: E0143 and E0156
- 5 Year Limited Warranty

RJ4700B Blue
RJ4700R Burgundy
RJ4700S Silver



Specifications*

Width

Inside Handgrips: 16.3" (414 mm)
Overall: 24.5" (622 mm)

Depth

At Minimum Seat Height: 25" (635 mm)
At Maximum Seat Height: 29" (737 mm)

Height

At Minimum Handle Height: 29" (737 mm)
At Maximum Handle Height: 38" (965 mm)

Seat

Height: 18" - 22" (457 mm - 546 mm)
Width: 14" (356 mm)
Depth: 12" (305 mm)
Casters Diameter: 6" (152 mm)
Maximum User Weight: 300 lb (136 kg)
Product Weight: 14 lb (6.4 kg)

*Specifications are subject to change without notice.



2 Rollators in 1!



GRAHAM-FIELD

colorSelect



www.grahamfield.com