

## Set n' Go Height Adjustable Rollator

The Lumex Set n' Go Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one. Designed with multiple height settings on the frame, the Set n' Go allows the seat height to be adjusted from a hemi 18" to a standard 22" height in 1" increments. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated.

- Easy height adjustment on the frame allows seat height to be adjusted from 18" to 22" in 1" increments
- Comfortable, wide, seamless padded seat with zippered pouch under seat for security of personal items
- Also includes black pouch for storage of larger items
- Ergonomic handgrips offer a comfortable grip
- Adjustable handle height with knobs to secure in place
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Removable padded backbar with easy release pull knobs
- 6" wheel for both indoor and outdoor use
- 300 lb maximum weight capacity
- Designed to meet the requirements of HCPCS: E0143 and E0156
- 5 Year Limited Warranty

**RJ4700B** Blue  
**RJ4700R** Burgundy  
**RJ4700S** Silver



### Specifications\*

#### Width

Inside Handgrips: ..... 16.3" (414 mm)  
Overall: ..... 24.5" (622 mm)

#### Depth

At Minimum Seat Height: ..... 25" (635 mm)  
At Maximum Seat Height: ..... 29" (737 mm)

#### Height

At Minimum Handle Height: ..... 29" (737 mm)  
At Maximum Handle Height: ..... 38" (965 mm)

#### Seat

Height: ..... 18" - 22" (457 mm - 546 mm)  
Width: ..... 14" (356 mm)  
Depth: ..... 12" (305 mm)  
Casters Diameter: ..... 6" (152 mm)  
Maximum User Weight: ..... 300 lb (136 kg)  
Product Weight: ..... 14 lb (6.4 kg)

\*Specifications are subject to change without notice.



**2 Rollators in 1!**



GRAHAM-FIELD

colorSelect



**www.grahamfield.com**