



## Aluminum Rollator, 7.5" Casters

Fold-Up and Removable Back Support, Flip-Up Seat, Removable Wheels, Loop Locks

Item #	Color	UOM
RTL728RD	Red 	1/cs

- Comes with new seamless padded seat. Zippered pouch under seat for added privacy and security of personal items. (Figure A)
- 7.5" Black non-marring removable casters.
- Ideal for indoor and outdoor use.
- Removable, hinged, padded backrest can be folded up or down.
- Easy-to-use deluxe loop locks. (Figure B)
- Brakes with serrated edges provide firm hold. (Figure C)
- Handles easily adjust in height with self-threading knob.
- Retail packaged.

### specifications

Handle Height	32.5" - 37"
Overall Dimensions	24" (W) x 29" (L)
Seat Dimensions	12" (W) x 23" (H) x 14" (D)
Weight	16 lbs.
Carton Shipping Weight	20 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

### accessories

9501S72606 Basket	1/ea
-------------------	------

## Aluminum Rollator, 7.5" Casters

Fold-up and Removable Back Support, Padded Seat, Loop Locks

Item #	Color	UOM
R728BL	Blue 	1/cs
R728RD	Red 	1/cs

- Comes with new seamless padded seat (Figure A). Zippered pouch under seat for added privacy and security of personal items. (Figure B)
- 7.5" Black non-marring tires are ideal for indoor and outdoor use.
- Removable, hinged, padded backrest can be folded up or down.
- Easy-to-use deluxe loop locks. (Figure C)
- Brakes with serrated edges provide firm hold. (Figure D)
- Ergonomic handles are easy to grip and are adjustable.

### specifications

Handle Height	33" - 38"
Unit Dimensions	24" (W) x 28" (L)
Seat Dimensions	14" (W) x 23" (H) x 12" (D)
Weight	14 lbs.
Carton Shipping Weight	17 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

### accessories

9501S72606 Basket	1/ea
-------------------	------