

RgoSling HighBack

Built-in ergonomy for comfortable hoisting.

For homecare and institutional environments

Molift RgoSling MediumBack and Molift Rgo Sling HighBack are all-round slings that fits most users and hoisting situations in homecare and institutional environments.

Molift RgoSling

RgoSling MediumBack provides support for the body and the RgoSling HighBack provides support for the head as well. It allows the user to have a slightly inclined position which is suitable for users with impaired trunk and head stability. It is especially suitable for hoisting to and from a lying position and hoisting from the floor.

Molift RgoSling in Net-version

The RgoSling MediumBack and RgoSling HighBack are also available in a soft net mesh polyester material. The material breathes well and dries very quickly. The sling is especially suitable for shower and bathing situations. The mesh version has the same shape and features as the padded versions, such as; sliding loops, application pocket in the back and straps for increased weight distribution in the leg parts.



RgoSling HighBack Net.

Molift RgoSling MediumBack

Size	Padded	Net
XXS		1720300
XS	1720110	1720310
S	1720120	1720320
M	1720130	1720330
L	1720140	1720340
XL	1720150	1720350
XXL	1720160	1720360

Molift RgoSling HighBack

Size	Padded	Net
XXS		1720400
XS	1720210	1720410
S	1720220	1720420
M	1720230	1720430
L	1720240	1720440
XL	1720250	1720450
XXL	1720260	1720460

Recommended User Weight

XXS	12–17 kg (26–37 lbs)
XS	17–25 kg (37–55 lbs)
S	25–50 kg (55–110 lbs)
M	45–95 kg (99–210 lbs)
L	90–160 kg (198–353 lbs)
XL	160–240 kg (353–529 lbs)
XXL	230–300 kg (507–660 lbs)

Material

Polyester,
PE -cellfoam

Washing instruction



Colour code

XXS	XS	
S	M	L
XL	XXL	

