







## **Blind Folding Cane**

Item #	UOM	
10352-1	1/cs	

- High quality four section aluminum construction (Figure A).
- Shaft covered with white and red reflective tape for night visibility.
- Cane length: 45.75".
- Comes standard with wrist strap and reinforced nylon tip.

## **Cane Seat**

Item #	UOM
10365	1/cs
RTL10365-ADJ	1/cs

- Provides a comfortable seat to rest on when open, and a sturdy support cane when closed. (Figure A)
- Manufactured with sturdy, extruded aluminum tubing.
- Tripod design with vinyl contoured tipped legs.
- 9" diameter seat.
- 33" tall when closed. (Figure B)
- Weighs 1.3 lbs.
- 250 lb. Weight Capacity.

specifications	
Cane Length	45.75"
Warranty	Limited Lifetime

specifications	
Seat To Floor Height	21"
Cane Height	34" - 38"
Seat Height	19" - 22"
Weight Limit	250 lbs.
Warranty	Limited Lifetime