



## Tall Extension Legs with Wheels, Combo Pack (Adds 4")

HCPCS:



FIGURE A



FIGURE B

- Allows for height adjustments from 36" to 43" (walker handgrip to floor) (Adds 4" to height of walker)
- 5" rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- For use with all Drive and most leading manufacturers walkers
- Includes rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces
- Weight Capacity 350 lbs

### PRODUCTS

#### Description

PRODUCTS	Description
10108WC	Ext.wheels Combo Pack,1 set/bx, RTL

### SPECIFICATIONS

#### Description

SPECIFICATIONS	Description
Weight Capacity	350 lbs