

The **Invacare® Storm Series® Arrow® Power Wheelchair** with TrueTrack technology and Invacare® MK6i™ Electronics has the flexibility and function to accommodate a variety of drivers with diverse needs.

Features

Adjustable rehab (ASBA) seat
 Captain's seat available in 18", 20" and 22" widths and 16"-20" depths
 Monroe® independent rear suspension
 Optional shock forks
 Transport Tie-Down options*
 Heavy-duty option available

Invacare Core Technologies

Invacare® Gearless Brushless GB™ Motors with TrueTrack technology
 Invacare® MK6i™ Electronics are the benchmark in innovative technology
 Motion Concepts power positioning systems can be interfaced with the Arrow power base. Contact Motion Concepts for more information at 888-433-6818

HCPCS Code: K0868, K0870, K0884, K0886



1. Invacare® Gearless Brushless GB™ Motor
2. Invacare® MK6i™ Electronics
3. Motion Concepts Powered Positioning System Compatibility



©2012

Invacare Corporation.
 All rights reserved.
 Trademarks are identified by the symbols ™ and ®.
 All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Monroe is a registered trademark of Tenneco automotive Inc. Specifications are subject to change without notification.
 Form No. 12-240
 120340



Model #	Overall Dimensions	Rehab Seat Dimensions	Rehab Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
3GAR	41.5" L	12"-16"x	17.25"-21.75"	14"	3"	31.9"	9°
3GAR-CG	26.6" W	12"-18"					

	Total Weight	Weight Capacity	Motor	Max Speed	Battery	Battery Range	Limited Warranty
3GAR	Up to 401 lb.	300 lb.	TrueTrack GB	6.7 mph	Group 24	Up To 29.2 Miles	5 Years Base Frame
3GAR-CG		400 lb.					3 Years Seat Frame 1 Year Electronics/Motor/Gearbox