



FlamingoTM

The Flamingo is a toilet/bathing chair. It is available in 4 sizes and is suitable for children and youngsters



The seat can be tilted and the back can be angled. The Flamingo can therefore assist special needs children with their seating posture. The Flamingo is easy to clean and facilitates personal hygiene. It is adjustable in height and it can be placed over the toilet or can be used with a pan.

As standard the Flamingo seat is delivered with a narrow open backrest. For more support for the trunk either swing-away side supports can be mounted or upholstered sides which surround and support the shoulders can be used.




Tilt in space and back recline

Tilts and reclines for the child's comfort and positioning



Adjustable foot supports

Adjust to ensure correct positioning and support



Height adjustable frame

Adjusts to fit the height of the carer



Foldable frame

Folds for storage or transportation

Sizes

88000x-NL

The Flamingo is available in four sizes



Cross bar

The cross bar is curved to fit around different commode designs



Curved

Rust-proof frame and locking casters




Aluminium frame
with lockable casters
for safe and stable
transfer

Frame & castors

The Flamingo can be used over a toilet, with a pan, directly on the toilet, as a bathing seat in the bath tub, as a shower chair or reclined over the bath for hair washing





Measurements

	Size 1	Size 2	Size 3	Size 4
Seat width (A)	25 cm (10")	28 cm (11")	32 cm (12½")	34 cm (13")
Seat depth (B)	24 cm (9½")	29 cm (11½")	34 cm (13¼")	40 cm (15½")
Back height (C)	36 cm (14")	38 cm (14¾")	41 cm (16")	47 cm (18½")
Lower leg length (D)	18-26 cm (7-10")	25-35 cm (9¾-13½")	30-40 cm (11¾-15½")	30-40 cm (11¾-15½")
Seat height (E)	50-60 cm (19½-23½")	50-60 cm (19½-23½")	50-60 cm (19½-23½")	50-60 cm (19½-23½")
Shoulder width (F)	30 cm (12")	32 cm (12½")	38 cm (15")	40 cm (15½")
Back width	19½ cm (7½")	22 cm (5½")	26 cm (10")	28 cm (11")
Width, frame (G)	44 cm (17")	46 cm (18")	50 cm (19¾")	54 cm (21")
Length, frame (H)	66 cm (26")	66 cm (26")	66 cm (26")	66 cm (26")
Height, frame w/seat (I)	86 cm (33½")	89 cm (35")	93 cm (36½")	95 cm (37")
Seat height, folded	24 cm (9½")	26 cm (10")	27 cm (10½")	27 cm (10½")
Seat length, folded	37 cm (14½")	40 cm (15½")	43 cm (17")	46 cm (18½")
Frame height, folded	20 cm (8")	20 cm (8")	20 cm (8")	20 cm (8")
Frame length, folded	74 cm (29")	74 cm (29")	74 cm (29")	74 cm (29")
Weight	11 kg (24 lb)	13 kg (29 lb)	14 kg (31 lb)	15 kg (33 lb)
Max. load/user weight	35 kg (77 lb)	50 kg (110 lb)	60 kg (132 lb)	70 kg (154 lb)
Angling	Seat:	-15°-+35°	Back:	-25°-90°
	Foot support bar:	90°-60°		

ACCESSORIES



Head support

9905145-00

For correct positioning and support

See the Accessories section in the back of the catalog for further information on options and measurements.



Push brace

880110

Provides secure point to maneuver the unit

Size	Measurements
1	W: 14"




Sides w/upholstery

88011x

Provide moderate torso control and positioning

Size	Measurements
1	W:3¾ x H:9"
2	W:3¾ x H:9½"
3	W:3¾ x H:10½"
4	W:3¾ x H:12½"



Side supports

83325-262

Swing away side supports for additional support and safety

See the Accessories section in the back of the catalog for further information on options and measurements.



Armrest/tray

88013x. Provides upper arm support and swings away for transfer

Size	Measurements
1	W:13½ x D:9"
2	W:13¾ x D:9"
3	W:16¾ x D:10¼"
4	W:18½ x D:10¼"



Armrest

880151

One size




Commode pan

88012x

For use inside or outside the bathroom

Size	Measurements
1+2	D:8¼ x H:7"
3+4	D:12¼ x H:7"



Upholstery, seat with wider sides

88016x. For comfort and reduced seat width - also making the standard hole narrower

Size	Measurements
2	Cushion W:1"
3	Cushion W:1"
4	Cushion W:1"



Fitting for toilet mounting

88014x. For mounting the Flamingo on an ordinary toilet

Size	Measurements
1	W:11½ x D:8"
2	W:12½ x D:10¼"
3	W:14½ x D:12"
4	W:15½ x D:14"



Net vest

89111-x40

For correct positioning and support

See the Accessories section in the back of the catalog for further information on options and measurements.



Calf support

880199

For added support and comfort

Size	Measurements
1	H:4 x L:42"




Splash guard, low

880195, 880196

Size	Measurements
1-4	Above seat, H:0"

ACCESSORIES



Suction cups

880150

For safe and stable mounting
in the bath tub

One size



Seat infill

88018x. For covering the
aperture in the seat

Size	Measurements
------	--------------

1	W:4 x D:15¼"
2	W:5 x D:7"
3	W:6 x D:9¼"
4	W:6 x D:11¼"



Wheels

880100-5.

Bigger wheels are available
for the Flamingo

4"
