

R82 Flamingo | Flamingo High-Low

CAD Retail Price List/ Order Form



Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____


| | |
|-----------------|-----------------|
| Bill To: | Ship To: |
| Address: | Address: |
| City/State/Zip: | |
| Phone: | City/State/Zip: |
| Email | Phone: |


Flamingo Standard Features:


- | | |
|--|--|
| <ul style="list-style-type: none"> Size 1, 2, 3, 4 4 White Locking Castors Aluminum Frame 100% Waterproof Upholstery for Seat and Back Upholstery Available in Blue | <ul style="list-style-type: none"> Height Adjustable and Foldable Frame Tilt in Space (+5° to -35°) Reclining Back Active Back (No Sides) Splash Guard - High Model |
|--|--|



Flamingo High-Low Standard Features:



- | | |
|---|---|
| <ul style="list-style-type: none"> Size 1, 2, 3, 4 4 White Locking Castors Foot Support Stainless Steel Frame 100% Waterproof Upholstery for Seat and Back Upholstery Available in Blue | <ul style="list-style-type: none"> Height Adjustable Frame Tilt in Space Reclining Back Active Back (No Sides) Splash Guard - High Model |
|---|---|



| | <u>Flamingo Seat and Back</u> | <u>Size</u> | <u>Part No.</u> |
|---|---|-------------|-----------------|
| (A)  | <input type="checkbox"/> Flamingo Seat and Back, Blue (A) | 1 | 880101-61 |
| | <input type="checkbox"/> Flamingo Seat and Back, Blue (A) | 2 | 880102-61 |
| | <input type="checkbox"/> Flamingo Seat and Back, Blue (A) | 3 | 880103-61 |
| | <input type="checkbox"/> Flamingo Seat and Back, Blue (A) | 4 | 880104-61 |



| | <u>Flamingo Frame</u> | <u>Size</u> | <u>Part No.</u> |
|--|--|-------------|-----------------|
| (A)  | <input type="checkbox"/> Flamingo Frame w/out Footrest (A) | 1 | 880201-NL |
| | <input type="checkbox"/> Flamingo Frame w/out Footrest (A) | 2 | 880202-NL |
| | <input type="checkbox"/> Flamingo Frame w/out Footrest (A) | 3 | 880203-NL |
| | <input type="checkbox"/> Flamingo Frame w/out Footrest (A) | 4 | 880204-NL |



| | <u>Flamingo High-Low Bath Frame</u> | <u>Size</u> | <u>Part No.</u> |
|---|--|-------------|-----------------|
| (A)  | <input type="checkbox"/> High-Low Bath Frame w/ Footrest (A) | 1 | 880411 |
| | <input type="checkbox"/> High-Low Bath Frame w/ Footrest (A) | 2 | 880412 |
| | <input type="checkbox"/> High-Low Bath Frame w/ Footrest (A) | 3 | 880413 |
| | <input type="checkbox"/> High-Low Bath Frame w/ Footrest (A) | 4 | 880414 |



| | <u>Back Option</u> | <u>Size</u> | <u>Part No.</u> |
|---|--|-------------|-----------------|
| (A)  (B)  | <input type="checkbox"/> Back Extension (H: 4" x L: 9") (A)* | 1, 2 | 880118-61 |
| | <input type="checkbox"/> Back Extension (H: 4" x L: 9" - 10.2") (A)* | 3, 4 | 880119-61 |
| | <input type="checkbox"/> Mounting Rod for Back Extension | 7.9" | 880120 |
| | <input type="checkbox"/> Back Reinforcement Kit (B) | 3, 4 | 880155 |
| (A)* Requires Mounting Rod for Back Extension | | | |


| | | <u>Pelvic Positioning</u> | <u>Size</u> | <u>Part No.</u> |
|---|--|--|---|-----------------|
| (A) |  | (B) |  | |
| | | | | |
| | | <input type="checkbox"/> Upholstery, Seat w/ Wide Sides, Blue (A)* | 2 | 880166-61 |
| | | <input type="checkbox"/> Upholstery, Seat w/ Wide Sides, Blue (A)* | 3 | 880167-61 |
| | | <input type="checkbox"/> Upholstery, Seat w/ Wide Sides, Blue (A)* | 4 | 880168-61 |
| | | <input type="checkbox"/> Seat Infill for Seat, Blue (B) | 1 | 880185-61 |
| | | <input type="checkbox"/> Seat Infill for Seat, Blue (B) | 2 | 880186-61 |
| | | <input type="checkbox"/> Seat Infill for Seat, Blue (B) | 3 | 880187-61 |
| | | <input type="checkbox"/> Seat Infill for Seat, Blue (B) | 4 | 880188-61 |
| | | <input type="checkbox"/> Seat Belt w/ Plastic Lock | One Size | 85100 |
| (A)* Does not include Splashguard. Reduces sz 2 & 3 width by 2", and sz 3 by 3.5" | | | | |

| | | <u>Foot & Leg Positioning</u> | <u>Size</u> | <u>Part No.</u> |
|-----|--|--|---|-----------------|
| (A) |  | (B) |  | |
| | | | | |
| | | <input type="checkbox"/> Flamingo Foot Support (A) | 1 | 880211 |
| | | <input type="checkbox"/> Flamingo Foot Support (A) | 2 | 880212 |
| | | <input type="checkbox"/> Flamingo Foot Support (A) | 3 | 880213 |
| | | <input type="checkbox"/> Flamingo Foot Support (A) | 4 | 880214 |
| | | <input type="checkbox"/> Calf Support (H: 3.9" x L: 41.7") (B) | One Size | 880199 |
| | | <input type="checkbox"/> Footrest Straps, Pair | | 85096 |

| | | <u>Trunk Positioning</u> | <u>Size</u> | <u>Part No.</u> |
|---|--|--|---|-----------------|
| (A) |  | (B) |  | |
| | | | | |
| | | <input type="checkbox"/> Bar, Straight, Stainless (A)** | 3¾" | 82569 |
| | | <input type="checkbox"/> Bar, Straight, Stainless (A)** | 6¾" | 82648 |
| | | <input type="checkbox"/> Cushion, Curved, Blue Each (B)** | A: 2½", B: 5¾" | 6075409 |
| | | <input type="checkbox"/> Cushion, Curved, Blue Each (B)** | A: 2½", B: 6½" | 6075408 |
| | | <input type="checkbox"/> Cushion, Curved, Blue Each (B)** | A: 3¾", B: 7½" | 6075407 |
| | | <input type="checkbox"/> Protection Pads, Each | | 89151 |
| | | <input type="checkbox"/> Swing-Away Bracket, Right, Each (C)** | A: 1¾", B: 4¾" | 83300-R |
| | | <input type="checkbox"/> Swing-Away Bracket, Left, Each (C)** | A: 1¾", B: 4¾" | 83300-L |
| | | <input type="checkbox"/> Swing-Away Bracket, Right, Offset, Each (D) | A: 2¼", B: 5" | 83500-R |
| | | <input type="checkbox"/> Swing-Away Bracket, Left, Offset, Each (D) | A: 2¼", B: 5" | 83500-L |
| | | <input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)* | 1 | 880111-61 |
| | | <input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)* | 2 | 880112-61 |
| | | <input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)* | 3 | 880113-61 |
| | | <input type="checkbox"/> Sides for Back w/ Upholstery, Blue (E)* | 4 | 880114-61 |
| | | <input type="checkbox"/> Vest, Water-Proof Mesh (F) | 0 | 89111-040 |
| | | <input type="checkbox"/> Vest, Water-Proof Mesh (F) | ¾ | 89111-07540 |
| | | <input type="checkbox"/> Vest, Water-Proof Mesh (F) | 1 | 89111-140 |
| | | <input type="checkbox"/> Vest, Water-Proof Mesh (F) | 2 | 89111-240 |
| | | <input type="checkbox"/> H-Harness w/ Plastic Lock | | 85105 |
| (A), (B), (C) ** Required to complete Swing-Away Lateral Assy | | | | |
| (E)* Cannot be combined w/ Swing-Away Lateral Supports | | | | |

| | | <u>Head Support</u> | <u>Size</u> | <u>Part No.</u> |
|-----|--|--|---|-----------------|
| (A) |  | (B) |  | |
| | | | | |
| | | <input type="checkbox"/> Head Support No. 1, Blue PU (A) | | 6071000 |
| | | <input type="checkbox"/> Head Support No. 2, Blue PU (B) | | 6072000 |
| | | <input type="checkbox"/> Head Support No. 3, Blue PU (C) | | 6073000 |
| | | <input type="checkbox"/> Head Support No. 4, Blue PU (D) | | 6074000 |
| | | <input type="checkbox"/> Turtle Bar System, Short (E) | 5.1" | 98619-24 |

| | | <u>Frame Options</u> | <u>Size</u> | <u>Part No.</u> |
|-----|--|--|---|-----------------|
| (A) |  | (B) |  | |
| | | | | |
| | | <input type="checkbox"/> Commode Pan w/ Holder for Flamingo (A) | 1 | 880121 |
| | | <input type="checkbox"/> Commode Pan w/ Holder for Flamingo (A) | 2 | 880122 |
| | | <input type="checkbox"/> Commode Pan w/ Holder for Flamingo High-low (A) | 1 | 880121-HL |
| | | <input type="checkbox"/> Commode Pan w/ Holder for Flamingo High-low (A) | 2 | 880122-HL |
| | | <input type="checkbox"/> Commode Pan w/ Holder for Flamingo & High-Low (A) | 3 | 880123 |
| | | <input type="checkbox"/> Commode Pan w/ Holder for Flamingo & High-Low (A) | 4 | 880124 |

Frame Options Continued

| | |
|--|----------|
| <input type="checkbox"/> Potty Ring w/ Holders (B)* | 1 |
| <input type="checkbox"/> Potty Ring w/ Holders (B)* | 2 |
| <input type="checkbox"/> Potty Ring w/ Holders (B)* | 3 |
| <input type="checkbox"/> Potty Ring w/ Holders (B)* | 4 |
| <input type="checkbox"/> Push Brace (C) | One Size |
| <input type="checkbox"/> Push Handle for Flamingo High-Low (D) | 1, 2 |
| <input type="checkbox"/> Push Handle for Flamingo High-Low (D) | 3, 4 |
| <input type="checkbox"/> Single Post Armrest (L: 9"), Set (E) | One Size |
| <input type="checkbox"/> Suction Cups for Seat (F) | One Size |

Part No.

| |
|------------|
| 880121-PR |
| 880122-PR |
| 880123-PR |
| 880124-PR |
| 880110 |
| 880110-HL |
| 880110-HLL |
| 880152 |
| 880150 |

(B)* Potty ring required for use over the commode

(E)* Armrest cannot be used together w/ Sides for Back - 88011X-X

Accessories

| | |
|--|---------|
| <input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A) | 1 |
| <input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A) | 2 |
| <input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A) | 3 |
| <input type="checkbox"/> Swing-Away Armrest/Tray, Blue (A) | 4 |
| <input type="checkbox"/> Splashguard, Low Model, Blue (B) | 1 |
| <input type="checkbox"/> Splashguard, Low Model, Blue (B) | 2, 3, 4 |

Part No.

| |
|-----------|
| 880131-61 |
| 880132-61 |
| 880133-61 |
| 880134-61 |
| 880195-61 |
| 880196-61 |

IMPORTANT: No returns or exchanges of hygiene products

Effective: 1/10/2020

Rev 5


Flamingo Configuration

- ☐ Flamingo/Flamingo High-low
- ☐ Head support
- ☐ Push handle/push brace
- ☐ Side supports
- ☐ Arm rest/tray
- ☐ Commode pan/potty ring
- ☐ Vests and belts



A (Back height) _____
 B (Arm height) _____
 C (Lower leg length) _____
 D (Seat depth) _____
 E (Chest width) _____
 F (Seat width) _____
 G (Height from seat to top of head) _____
 H (Shoulder width) _____
 Age _____
 Weight _____

Measurements



| | Size 1 Inches | Size 2 Inches | Size 3 Inches | Size 4 Inches |
|--------------------------|------------------|------------------|------------------|------------------|
| Seat width (A) | 9¼" | 11" | 12½" | 13¾" |
| Seat depth (B) | 9¾" | 11¼" | 13¾" | 15½" |
| Back height (C) | 13¾" | 14¾" | 16½" | 18¼" |
| Lower leg length (D) | 18" - 26" | 25" - 35" | 30" - 40" | 30" - 40" |
| - | 7" - 10" | 9¾" - 13½" | 11¾" - 15½" | 11¾" - 15½" |
| Seat height (E) | 50" - 60" | 50" - 60" | 50" - 60" | 50" - 60" |
| - | 19½" - 23½" | 19½" - 23½" | 19½" - 23½" | 19½" - 23½" |
| Back width (F) | 7¾" | 8½" | 10¼" | 11½" |
| Width, frame (G) | 17" | 18" | 19¼" | 21" |
| Length, frame (H) | 26" | 26" | 26" | 26" |
| Height, frame w/seat (I) | 33½" | 35" | 36½" | 37" |
| Seat height, folded | 9½" | 10" | 10½" | 10½" |
| Seat length, folded | 14½" | 15½" | 17" | 18½" |
| Frame height, folded | 8" | 8" | 8" | 8" |
| Frame length, folded | 29" | 29" | 29P | 29" |
| Weight | 24 lbs. | 29 lbs. | 31 lbs. | 33 lbs. |
| Max. Load/user weight | 77 lbs. | 110 lbs. | 132 lbs. | 154 lbs. |
| Angling | | | | |
| Seat | | -15°-+35° | Back | -90°-30° |
| Foot support bar | | | 90°-60° | |

Flamingo High-Low Configuration

- ☐ Flamingo/Flamingo High-low
- ☐ Head support
- ☐ Push handle/push brace
- ☐ Side supports
- ☐ Arm rest/tray
- ☐ Commode pan/potty ring
- ☐ Vests and belts



A (Back height) _____

B (Arm height) _____

C (Lower leg length) _____

D (Seat depth) _____

E (Chest width) _____

F (Seat width) _____


G (Height from seat to top of head) _____

H (Shoulder width) _____

Age _____

Weight _____

Measurements



| | Size 1 Inches | Size 2 Inches | Size 3 Inches | Size 4 Inches |
|-----------------------------|------------------|------------------|------------------|------------------|
| Seat width (A) | 9¼" | 11" | 12½" | 13¾" |
| Seat depth (B) | 9¾" | 11¼" | 13¾" | 15½" |
| Back height (C) | 13¾" | 14¾" | 16½" | 18¾" |
| Lower leg length (D) | 7¼" - 9¾" | 9" - 13¼" | 11¾" - 16½" | 11¾" - 16½" |
| Seat height above floor (E) | 20¾" - 28¾" | 20¾" - 28¾" | 20¾" - 28¾" | 20¾" - 28¾" |
| Back width (F) | 7¼" | 8½" | 10¼" | 11" |
| Frame width (G) | 19½" | 19½" | 22¼" | 22¼" |
| Frame length (H) | 25" | 25" | 29½" | 29½" |
| Total height (I) | 33½" - 42" | 35¾" - 43¼" | 36¾" - 45¾" | 39¾" - 46¾" |
| Seat, total width | 13¾" | 14¾" | 16¾" | 18¾" |
| Seat, total depth | 10¼" | 12¼" | 14" | 16½" |
| Seat, total height | 15½" | 16¾" | 18¾" | 19¾" |
| Back angling | -90°-+30° | -90°-+30° | -90°-+30° | -90°-+30° |
| Tilt in space | -5°-+20° * | -5°-+20° | -5°-+20° | -5°-+20° |
| Foot support angle | 90°-75° | 90°-75° | 90°-75° | 90°-75° |
| Weight, total | 31 lbs. | 32 lbs. | 36½ lbs. | 37½ lbs. |
| Weight, frame | 25½ lbs. | 25½ lbs. | 27½ lbs. | 27½ lbs. |
| Weight, seat | 5½ lbs. | 6½ lbs. | 9 lbs. | 10 lbs. |
| Max. load, total | 77 lbs. | 110 lbs. | 132 lbs. | 154 lbs. |
| Max. load, seat | 77 lbs. | 110 lbs. | 132 lbs. | 154 lbs. |
| Max. load, operating force | 77 lbs. | 77 lbs. | 121 lbs. | 121 lbs. |

* The -5° tilt is obtainable at 55 cm (19½") seat height above floor