

Dolphin End-Loading Bath Chair

Bathing in seated comfort

Made in Canada



Equipment Solutions for Patient Bathing & Lifting



The unique, fully integrated scale offers a simple and reliable means of weighing patients / residents.



An optional retractable / removable leg support is also available.

End-Loading Bath Chair

The Dolphin End-Loading Bath Chair's ergonomic push handle makes it easy to manoeuvre. The end-loading design makes room to room transfers safer and easier than other designs. With an extensive range of operating heights, the chair can be used with a variety of tubs, including the Dolphin and Porpoise bath systems. Other systems must have an inside slope for the backrest of 12 degrees and be clear underneath to floor to accommodate the chair legs.



Model #DBC-208IS

Features

- Wide range of height adjustments from 20.75" to 42.5" (52.7 cm to 110.5 cm)
- Stainless Steel and durable epoxy Silver Vein coated surfaces
- Retractable wrap-around safety arms for patient security and easy entry
- Adjustable safety belt included
- IBS (Interchangeable Battery System - see below)

Seat and Backrest

- Designed to be comfortable, durable and easy to clean
- Material is: tear resistant, quick drying, fluid proof, low shear, stain resistant, anti-bacterial, anti-microbial, fungus resistant, anti-static, latex-free, and radiolucent



IBS Power Pack

Used on all Care Source lifting devices, retrofit kits, chairs and stretchers, this fully interchangeable system simplifies use and creates continuity of operation.

- Modular components for at source servicing
- Wall-mounted exchange pack charging station
- Durable hand control
- Emergency stop button and lowering switch
- CSA and UL approved charger

For more information, visit www.caresourcetek.com



CANADA: Care Source Inc., 2026 Headon Road, Burlington, ON L7M 4G3
U.S.A.: Care Source Technologies Inc., 485 Five Points Road, Alburts, PA 18011
 T: 905.635.1169 • F: 289.337.5337
www.caresourcetek.com info@caresourcetek.com