

Dolphin Side-Loading Bath Chair

Bathing in seated comfort



Equipment Solutions for Patient Bathing & Lifting



The unique, fully integrated scale offers a simple and reliable means of weighing patients / residents.



An optional retractable / removable leg support is also available.

Side-Loading Bath Chair

The Dolphin Side-Loading Bath Chair's ergonomic push handle makes it easy to manoeuvre. With an extensive range of operating heights, the versatility of this side-loading chair allows it to work with many different bathing systems, including competing brands. The arms, backrest and seat easily change position to allow tub access from either side.

Features

- Wide range of height adjustments from 20.75" to 42.5" (52.7 cm to 110.5 cm)
- Stainless Steel and durable epoxy SilverVein coated surfaces
- Easily interchangeable backrest, arm and seat allow lift from either side
- Adjustable safety belt included
- IBS (Interchangeable Battery System - see below)

Seat and Backrest

- Designed to be comfortable, durable and easy to clean
- Material is: tear resistant, quick drying, fluid proof, low shear, stain resistant, anti-bacterial, anti-microbial, fungus resistant, anti-static, latex-free, and radiolucent



Model #DBC-208IS



IBS Power Pack

Used on all Care Source lifting devices, retrofit kits, chairs and stretchers, this fully interchangeable system simplifies use and creates continuity of operation.

- Modular components for at source servicing
- Wall-mounted exchange pack charging station
- Durable hand control
- Emergency stop button and lowering switch
- CSA and UL approved charger

For more information, visit www.caresourcetek.com



CANADA: Care Source Inc., 2026 Headon Road, Burlington, ON L7M 4G3 T: 905.635.1169 • F: 289.337.5337
U.S.A.: Care Source Technologies Inc., 485 Five Points Road, Alburtis, PA 18011 T: 610.682.1081 • F: 289.337.5337

www.caresourcetek.com info@caresourcetek.com