

# SPRITE SERIES.

The Sprite series is the most economical choice in the Evolution product line without compromising any quality or reliability.

Available in three sizes including two models with 40mm wide 6" wheels, suitable for indoor and outdoor use. The Sprite Grande comes with fork braking system, 8" wheels and is a lightweight full size walker designed for full indoor and outdoor use.

## FEATURES.

**Innovative Design** Light weight & easy to transport  
**Seat** Solid plastic base with removable fitted waterproof cushion  
**Wire Basket** Removable  
**Back Support** Solid padded  
**Under Seat Bag** (Sprite Regular Only) Provides extra storage  
**Colour** Vibrant Red, Sapphire Blue, Silver, Black



## SPECIFICATIONS.

Mini



11 lbs  
17" Seat Width  
18" Seat Height  
24" Width x 22.5" Length  
30" - 36" Handle Height  
250 lbs Weight Capacity

Regular



12 lbs  
17" Seat Width  
21" Seat Height  
24" Width x 23.5" Length  
33" - 38" Handle Height  
250 lbs Weight Capacity

Grande



13 lbs  
17" Seat Width  
21" Seat Height  
24" Width x 25.5" Length  
33" - 38" Handle Height  
250 lbs Weight Capacity

## OPTIONAL ACCESSORIES.

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes (Sprite Mini & Regular)
- Oxygen Tank Holder
- Tray - Seat Mounted
- Tote Bag (fits inside Mini & Grande basket)



**Evolution**  
TECHNOLOGIES\*

EVOLUTION TECHNOLOGIES INC.  
2530 Davies Avenue,  
Port Coquitlam, BC V3C2J9

tf +1 800.556.2558  
info@evolution.ca  
www.evolution.ca

EVOLUTION TECHNOLOGIES USA INC.  
2657 Delta Ring Road  
Ferndale WA 98248

tf +1 800.839.2558  
info@evolutionwalker.com  
www.evolutionwalker.com