

FORTRESS

Nothing is Built like a Fortress

LUXURY IN MOTION

S700

S700

PERFORMANCE LUXURY SCOOTER

The Fortress S700, with superior performance, features the perfect companion for all your adventures.

There's no need to worry about distance either as the large battery capacity will power an extended range and a lively top speed of **15 km/h**.

The high performance design means varied terrains are tackled with ease, and the large tires and advanced 4-wheel suspension provides a smooth, comfortable ride each and every time.





Max Speed: **15 km/h**



Colour: **Carbon Metallic Grey**



Max. Battery Capacity: **75 Amps (GP24)**



Wheel Size: **33 cm (13 inches)**



Max. User Weight: **159 kg. (350 lbs.)**

A close-up photograph of a scooter's front end. The main headlight is illuminated, casting a bright white glow. Below it, a smaller, rectangular turn signal is also illuminated, casting a bright yellow glow. The scooter's body is dark, and the background is black.

DELUXE LED LIGHT PACKAGE

**Powerful all-around
LED lighting. With
lifetime performance,
changing bulbs is a
thing of the past.**

INNOVATION IN DESIGN

With advanced LED lighting technology, the S700 scooter is 400 times more energy efficient than standard bulbs, meaning more battery power and range every trip.


**Safety is standard! S700
features high visibility LED turn
signals, 4 way flashers and
automatic brake lights.**



Innovative active anti-tip wheels
move up as you ride over obstacles, and
keep you safe on slopes.

LUXURY SEATING

**The multi-adjustable luxury seating
is designed to make every journey a
comfortable one.**



COMFORT ON EVERY TRIP

The S700 advanced seating solution takes comfort to the highest level. From the seat height, seat depth and back recline adjustments, to the flip-up armrests including width, angle and depth adjustments – every feature has been designed for comfort – for you.



→ Infinite seat depth adjustment



→ Quick backrest angle adjustment



→ Seat rotation enables comfortable transfers

CONTROLS DESIGNED FOR YOU

Finding your perfect driving position is a breeze with the infinitely adjustable tiller. Clever design means the headlight automatically tilts as the tiller is adjusted.

Soft touch controls featuring easy to use speed adjustment, LED battery gauge and a powerful horn.



Unique forward/reverse lever design means natural driving using finger or thumb control. The large ergonomically designed delta handlebar makes every journey a comfortable one.

Safety and style are standard. Every S700 scooter comes with two carbon effect mirrors.



Deluxe extra large basket comes standard to ensure plenty of room for your belongings.

S700

TECHNICAL DATA

	<i>S700</i>
Max Speed	15 km/h (9.3 mph)
Max. User Weight	159 kg. (350 lbs.)
Max Battery Capacity	75 Amps (GP24)
Charger	8 A off-board
Max. Potential Range*	55 km (34 miles)
Deluxe LED Light Package	Yes
Easy Adjust Tiller & Throttle Control	Yes
Suspension	Deluxe Front & Rear Suspension (Front Independent suspension)
Tires	13" pneumatic
Overall Width / Length	66 cm / 149 cm (26" / 59")
Seat Width x Depth	48 cm x 48 cm (19" x 19")
Seat Height Range	45 cm - 50 cm (18" - 20")
Minimum Turning Distance	2.2 m (7 feet)
Colour	Carbon Metallic Grey
Scooter Weight w/ Batteries 75 Amp	320 lbs. inc. GR24

* Max potential range based on industry benchmark range figures for similar products and battery sizes

All over Canada, thousands of people have made their lives easier with a FORTRESS mobility scooter by Sunrise Medical. For over 30 years, FORTRESS has designed and built mobility scooters to the highest standards, giving the comfort and safety every customer needs. For further information on product options please check the order form on our website. You can view this document also in a larger font on a PDF file at www.SunriseMedical.ca



Sunrise Medical Canada Inc.

237 Romina Drive, Unit 3
Concord, ON

Tel.: 1.800.263.3390

Fax: 1.800.561.5834

E-mail: mkt.canada@sunmed.com

www.SunriseMedical.ca