



Note: cover and foam insert required for use.

ROHO® Hybrid Elite SR™ Single Compartment Cushion

This innovative cushion combines the skin and soft tissue protection from the Shape Fitting Technology® of an air-filled ROHO DRY FLOATATION® Cushion with stability from a contoured foam base. The shape-fitting air insert provides critical pressure redistribution around the pelvis area.

With the available Smart Check, you can inflate your cushion to the proper level – and check it when ever you would like with the push of a button, minimizing tissue deformation and facilitating blood flow.



Personal inflation monitoring



Skin/soft tissue protection



Pressure injury prevention and treatment



Transfer stability

Dimensions

Height: 4.25 in. / 11 cm
Width: 13.75 in. to 23.75 in. / 35 cm to 60.5 cm
Weight: Varies by size, approximately 3 lbs. / 1.4 kg

Weight limit

500 lbs. / 225 kg The cushion must be properly sized to the individual.

Limited warranty

Cushion: 24 months
Cover: 6 months



Our evidence shows the Hybrid Elite SR is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require accommodation of pelvic asymmetry up to 2 in. (5 cm)
- Require a more stable environment for transfers
- Require pelvic and lower extremity positioning
- Require increased sitting stability