

Ordering Information

Part Number Description

Cefar Rehab X2

2942	Rehab X2, Complete kit including hand switch (EN-ES-FR-IT-SE-DE-PT-FN-TU)
------	---

Standard Accessories

1240	Set of Lead Wires (2 pcs)
2948	Cefar Rehab X2 User Manual (EN-ES-FR-IT-SE)
1559	Cefar Bag (Black), 2ch
1495	Cefar Neck Strap
1551	Batteries, Alk 1.5V AA (2 pcs)
1845	Hand Switch

Optional Accessories

1570	Fast Charger 12V
1575	Rechargeable Batteries 1800 mAh

For Electrodes, please see pages 17-18

Technical Specifications

Number of Channels:	2
Number of Programs:	27 preset, 3 programmable
Constant Current:	Up to a resistance of Max 1000 Ohm/channel
Maximum Amplitude:	99.5 mA
Intensity Increments:	Minimum 0.5 mA
Max Pulse Duration/Width:	400 μ seconds
Frequency:	1-120 Hz
Power Source:	2x1.5V AA
Dimensions*:	120 x 30 x 50 mm
Weight:	180g, including 2 AA batteries
Electrical Safety Class:	Class II, Type BF
Safety Tests:	IEC 60601-1, IEC 60101-1-2

* Dimensions are listed as L x W x H

Waveforms



- Symmetrical Biphasic 100% compensated

Cefar® Rehab X2



Cefar Rehab X2

The Cefar Rehab X2 offers 30 programs covering a large number of indications.

Features

- TENS: Cefar Flow TENS, mixed frequency (Han-stimulation), several pain relief programs
- NMES: Rehab programs for neurological and orthopaedic rehabilitation
- Alternating stimulation, strength programs with active rest, warm up, recovery and massage programs

Technology

- Cefar Easy Touch™: Automatic Intensity Adjustment - At the press of a button, the intensity is automatically set and adjusted according to the results of the Auto test analysis Cefar Easy Touch makes treatment easier and ensures stimulation takes place at an optimal and effective level
- Programmable: Cefar Rehab X2 allows the therapist to create and save three personal programs

- Cefar Flow Tens: A preset program featuring a comfortable stimulation method. Use both channels at the same time to combine pain relief with a pleasant massaging sensation. Other advantages: backlit screen, pause function, program lock, timer, key lock, self check and compliance
- Two year warranty

Advantages

- Cefar Clinical Guide (CCG™): Comprehensive information about the most common indications can be found in the user manual. Find guidelines and suggested programs with photos of electrode placements for the most common indications
- Hand Switch: The hand switch allows manual adjustment of the contraction time in intermittent stimulation programs, thereby ensuring effective and personalized treatment

