



## Nitro HD, Black

Item #	Product Description	UOM
10266HD-BK	Nitro HD Rollator, Black	1/ea

- Attractive, Euro-style design
- Brake cable inside frame for added safety (Figure A)
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort (Figure B)
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage (Figure C)
- Under seat locking mechanism holds the Nitro closed when frame is folded (Figure D)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

### specifications

Handle Height	36"-41"
Unit Dimensions	28" (L) x 28.25" (W)
Seat Dimensions	23" (W) x 10" (D)
Seat Height	23.5"
Weight	20.5 lbs.
Weight Limit	450 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake Handle & Housing)	5 Years

### accessory

892	Cane Holder	1/ea
10266-CH	Nitro Cup Holder	1/ea