wenzelite





Bariatric Anterior Safety Roller

ltem #	Color	Product Width Inside Handle Grip	Product Dimensions (open at base)	Product Weight	Product Weight Capacity	UOM
CE 1000 XL	Silver Vein 📕	17.5" - 22.5"	22.5" - 27.5" (W) x 25" (D)	19 lbs.	500 lbs.	1/bx
CE OBESE XL	Black	19" - 25"	24" - 30" (W) x 25" (D)	29 lbs.	1000 lbs.	1/bx

- Welded steel frame.
- Height adjustable in 1" increments.
- Width adjustable. (Figure A)
- Adjustable-brake spring tension.
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture.
- Skid-resistant, non-marring, easy-rolling neoprene wheels. Ideal for indoor and outdoor use.
- Folds easily and stands on its own in the folded position. (Figure B)

specifications	
Handlebar Height	29" - 36"
Carton Shipping Weight (CE 1000 XL)	24 lbs.
Carton Shipping Weight (CE OBESE XL)	35 lbs.
Warranty	Limited Lifetime