

manateeTM

The Manatee is constructed to provide both caregivers and users a wonderful bathing experience



The bathing chair is designed to place the user in exactly the right sitting position. The depth of the seat can be adjusted, allowing the bathing chair to "grow" with the child and provides many years of use. The adjustments and accessories are logically positioned, affording the carer better working conditions.

The Manatee is a bathing chair suitable for children from the age of 1 year. Not all children are alike - that is why it is important to use the measuring chart to find the correct size for your child.



Depth adjustable seat

Provides depth adjustment for growth



Adjustable back and leg rests

Allow adjustments of the back and of the lower leg support with a simple pull of a strap or push of a button



Head cushion/lateral support

Provides support and comfort for bathing children. Cushions can be used for head or torso support



Adjustable waist strap

Supports the child safely in the Manatee

Sizes


8820xx

The Manatee is available in four sizes




Cover colors

The Manatee is available in two different colors. Size 0 is only available with blue cover



Red




Blue

Height adjustment and variations



The Manatee can be lowered to within 5" of bottom of the bath. It is lightweight and it folds flat for storage. Also the seat can be placed directly into an ordinary bath tub




Measurements

	Size 0	Size 1	Size 2	Size 3
Seat width (A)	33 cm (12¾")	33 cm (12¾")	33 cm (12¾")	40 cm (15½")
Seat depth (B)	22 & 25 cm (8½-9¾")	30-37 cm (11¾-14½")	39-46 cm (15¼-18")	40-47 cm (15½-18¼")
Back height (C)	43 cm (16¾")	43 cm (16¾")	58 cm (22½")	64 cm (25")
Lower leg length (D)	28 cm (11")	28 cm (11")	35 cm (13¾")	35 cm (13¾")
Seat height (E)	14-43 cm (5½-16¾")	14-43 cm (5½-16¾")	14-43 cm (5½-16¾")	14-43 cm (5½-16¾")
Total width (G)	42 cm (16½")	42 cm (16½")	42 cm (16½")	49 cm (19")
Length folded	103-110 cm (40-43")	110-117 cm (43-45¾")	136-143 cm (53-55¾")	144-151 cm (56-58¾")
Height folded	13 cm (5")	13 cm (5")	13 cm (5")	13 cm (5")
Weight	5½ kg (12 lb)	5½ kg (12 lb)	6 kg (2¼ lb)	7 kg (15½")
Max. load	50 kg (110 lb)	50 kg (110 lb)	75 kg (165 lb)	75 kg (165 lb)
	Height adjustable frame		Basic frame	
Frame width	59 cm (23")		43½ cm (17")	
Frame length	98 cm (38¼")		99-111 cm (38½-43¼")	
Frame height min, mid, max	38-44-48 cm (14¾-17-18¾")		27 (10½")	
Frame height folded	21 cm (8")		-	
Frame weight	8 kg (17,5 lb)		2 kg (4,5 lb)	
Max.load on frame	75 kg (165 lb)		75 kg (165 lb)	

ACCESSORIES




Torso strap
883025 or 883075
For correct position and support



Pommel
883045 or 883095
For better positioning of the legs and prevents leg abduction



Leg supports
883055 or 883105
Provide precise control and positioning for the legs



Pads for head/lateral support

Size
0+1+2 88301x
3 88306x



Frame - Basic version
882101
For a better working position for the helper. To be placed in a bathtub



Frame - height adjustable
882100
For a better working position for the helper. To be placed in the shower



Bath top slider
For a better working position for the helper. To be placed in a bathtub



Net
88210x
For the height adjustable frame



Wall mounting
882110
For the height adjustable frame